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Deguform

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso

Dent

DeguDent

Deguform®

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Deguform®

1. Description

Deguform is a high-quality addition-curing silicone duplicating compound for producing secondary models in dental technology.

The two-component silicone is suitable for all modern duplicating techniques and convinces by:

- excellent reproduction of detail
- form-stable accurate reproduction
- very good tear resistance
- flow-viscosity consistency
- firm and smooth model surfaces
- simple and safe handling by 1:1 dosage
- time saving opposite gel duplication

2. Processing notes

2.1 Deguform is particularly suitable for duplicating in duplicating flasks. For this purpose prepare the Deguform duplicating flask – consisting of flask base and flask housing, individually adjustable palate and foam material sleeve.

2.2 Do **not** water the master model but duplicate dry and clean.

Note:

Silicone duplicating compounds are sensitive opposite catalyst poisons. These can be: e.g. insulating means, cyanoacrylate, plastics. Curing can thereby be delayed or prevented.

Lubricating the duplicating flask can also result in a non-hardening smear layer.

2.3 Place the foam liner around the master model to save on duplicating material. Coated side upwards.

Place the such prepared model on the flask base. Join palate and flask housing and place on the flask base.

Then adjust the palate to the individually required duplication height; it is held in this position for constructional reasons.

2.4 Mixing Deguform

Deguform is in principle mixed in the ratio of 1:1 (Component A, white: Component B, blue).

Deviations can poss. affect curing as well as the physical properties.

It is possible to mix Deguform using different methods:

- **by hand**
- **with suitable dosing devices**
- **with the Multivac® stirring units**


2.4.1 Mixing by hand

Fill Deguform in a ratio of 1:1 into a dosing and mixing vessel and mix thoroughly until a homogeneous colour results. Then immediately pour the compound in a thin stream into the duplicating flask.

2.4.2 Use of dosing devices

The minimum filling level must be observed during the use of 6 kg cans in connection with mixing and dosing units. Non-observance during processing results in air being fed into the system, a homogeneous mixture is no longer ensured. In such a case wait after connecting a new can for at least 45 minutes prior to the next withdrawal.

The following illustration shows the minimum silicone level in the can:

2.4.3 Mixing in Multivac® units

Have a separate mixing vessel ready for mixing silicone duplicating compounds. You achieve the best results with a new, unscratched vessel.

Fill Deguform in a ratio of 1:1 into the Multivac mixing unit, briefly stir thoroughly and mix in the Multivac according to the data in the following table. Then immediately pour the compound in a thin stream into the duplicating flask. The duplicating silicone can be easily be taken out of the mixing bowl and off the stirring blade after curing.

You must follow the recommended mixing parameters!

These have been determined so as to protect your equipment and to obtain optimal mixing results. Exceeding the maximum recommended **filling quantity** can in particular cause cementing of the filter and in extreme cases a drawing-in into the vacuum pump.

Device type	Bowl size	Max. filling quantity	Vacuum	Stirring speed	Stirring time
Multivac S, 4	Small	100 g	On	Presetting	60 s
Multivac S, 4	Large	400 g	On	Presetting	60 s
Multivac Eco	Small	100 g	On	Presetting	75 s
Multivac Eco	Large	400 g	On	Presetting	75 s
Multivac compact	Small	100 g	90 %	300 1/min	60 s
Multivac compact	Large	400 g	90 %	300 1/min	60 s

We recommend starting the stirring procedure only after achieving the vacuum. Longer evacuating does not affect the mixing result.

The silicone should always be filled to the extent that the three palate openings of the duplicating flask are filled and serve as an additional retention.

- 2.5 The silicone is cured after at least 30 minutes (at room temperature). Turn the duplicating flask over. Take flask base off. Remove foam liner. Remove the master model from the silicone form (if necessary with compressed air).

NOTE!

Removing the master model causes a deformation of the silicone form. This is neutralised after some minutes by the inherent resilience. We therefore recommend to only begin producing the investment model after a 15-minute waiting period.

- 2.6 The silicone form is highly elastic and keeps the precise form permanently.

Moisten the silicone form with Waxit® and blow dry before filling with the investment compound. The investment compound flow is thereby facilitated on the silicone surface. Make sure that the silicone surface is dry and completely set.

All kinds of dental gypsum and model plastics are suitable for filling additionally to investment compound.

- 2.7 Deguform can be disinfected both chemically and thermally.

3. Technical data

Component A	white
Component B	blue
Final colour	homogeneous light blue
Mixing ratio	1:1
Mixing time (manual)	approx. 45 s
Processing time *	approx. 5 min
Demoulding capacity*	approx. 30 min
Final hardness (shore A)	14 – 16

* Room temperature of approx. 23 °C. The time data can change at extreme room temperatures.

4. Storage conditions

Deguform should be stored at room temperature; only then can a minimum durability of 2 years starting from date of manufacture be ensured.

5. Product range

Deguform silicone duplicating compound, Component A and component B	REF
each in a bottle with 860 ml (1,000 g) contents	2511 0002
each in a can with 5,170 ml (6,000 g) contents	2511 0012
each in a can with 18.1l (21 kg) contents	2511 0042

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DeguDent GmbH