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# Optivest SR<sup>®</sup>

Instructions for use  
Istruzioni per l'uso  
Instrucciones para el uso  
Mode d'emploi

# Optivest SR®

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Ihr DeguDent Team

# Optivest SR®

## The multi-talent for lost-wax castings and non-precious alloys

**Optivest SR** is a gypsum-free, phosphate-bound and **low-dust** investment compound for duplicating casts and casting moulds.

A constantly optimized distribution of particles of different granularity yields extremely smooth casting surfaces.

Stable cast surfaces and exact reproduction of details ensure an excellent fit.

## Quality throughout

Within the framework of a certified quality management system, each batch is checked for correct product properties during the manufacturing process. The final product, too, is subjected to strict quality checks before being released for delivery to the customer.

This gives the dental technician a reliable investment compound that will yield top results, again and again!

## Possible applications

**Optivest SR** is pre-heated conventionally using holding levels. Optivest SR can be used with agar-agar and with silicone duplicating moulds.

### Manufacturer:

DeguDent GmbH

Rodenbacher Chaussee 4

63457 Hanau-Wolfgang

Germany

Phone +49 (0)61815950

Fax +49 (0)6181595858

### © December 2011 by DeguDent GmbH

We reserve the right to make changes and correct errors. This documentation may not be duplicated without our permission.

# Contents:

- 1 Safety information . . . . . 14
  - 1.1 Information for the processing of mixing liquids for dental investment materials . . . . . 14
- 2 Technical specifications . . . . . 14
- 3 Important information . . . . . 15
  - 3.1 Investment storage . . . . . 15
- 4 Fabrication of duplicate casts . . . . . 16
  - 4.1 Preparation . . . . . 16
  - 4.2 Mixing the investment compound . . . . . 16
  - 4.3 Pouring silicone duplicating moulds . . . . . 16
  - 4.4 Pouring agar-agar duplicating moulds . . . . . 16
  - 4.5 Preparation and drying . . . . . 17
  - 4.6 Hardening/Preparation . . . . . 17
- 5 Creating the casting mould . . . . . 18
  - 5.1 Stress relief . . . . . 18
  - 5.2 Protective investment coating . . . . . 18
  - 5.3 Types of flask . . . . . 18
  - 5.4 Mixing the investment compound . . . . . 18
  - 5.5 Investing . . . . . 18
  - 5.6 Hardening . . . . . 18
  - 5.7 Wax elimination . . . . . 18
  - 5.8 Pre-heating . . . . . 19
  - 5.9 Casting . . . . . 19
  - 5.10 Divesting . . . . . 19
- 6 Expansion control . . . . . 19

## 1 Safety information

Mixing with other liquids or materials during processing health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

### 1.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

## 2 Technical specifications

Mixing ratio (powder : liquid)	100 g : 15 ml
Mixing time (under vacuum)	60 s
Processing time	3 – 4 min
Setting time (Vicat)	4 – 6 min
Total linear expansion = setting expansion + thermal expansion*	0.7 – 1.8 %
Compressive strength**	13–20 MPa

\* Liquid concentration between 0 % and 100 %

\*\* Liquid concentration: 100 %

### 3 Important information

Be sure to observe information concerning safety prior to use under sections 1 and 1.1.

**Optivest SR** must not come into contact with plasters or plaster-containing materials since this may lead to sulphur damage of heat melting-alloys. Therefore special mixing bowls and spatulas must be used for mixing **Optivest SR**.

#### 3.1 Investment storage

Most of you are probably aware of the problems which can occur when investment is subjected to variations in environment temperature (e.g. summer and winter) and the effect this can have on your castings.

That is why it is very important to find the right place to store your supply of investment. The investment should be stored in a cool and dry place, between +5 °C to +30 °C at 30 % to 90 % relative humidity. If you have a place in your laboratory where these conditions can be found throughout the year, then that's the place to store your Investment!

If it is not possible for you to store your investment at this temperature, then you should at least make sure that the work bench where investing is done fits these requirements (use a temperature control cupboard or similar).

*Please note the symbols on product labels:*

-  REF Product code
-  LOT Batch or lot number
-  Use before
-  Consult instructions for use

*Avoid any nasty surprises by observing the following:*

- Regularly check the temperature and humidity in the embedding region - application temperature: +20 °C to +23 °C at 30 % to 90 % relative humidity.
- Avoid using any investment that has just been brought out of an especially warm or cold place. (e.g. storage rooms, silos, cupboards etc.)
- Fill up stocks at your workbench before your supply runs out.
- Give the investment time to adjust to the environment!
- Avoid using investment that has just been delivered to the laboratory. Again, give the investment time to acquire the right temperature.
- If space allows, stock up on investment (especially LIQUID) before the cold winter months arrive.
- This way, you forego the risk of exposing the liquid to frost during delivery.

## 4 Fabrication of duplicate casts

### 4.1 Preparation

Our motto is: "Clean tools – clean results!"

**Practical hint:** Use only clean mixing bowls that have been slightly moistened. Residual investment compound and scratches on the walls of the mixing bowl can trap large amounts of mixing liquid, adversely affecting not only the working time but also the fit.

### 4.2 Mixing the investment powder

**Mixing ratio (powder : liquid) = 100 g : 15 ml**

Rapidly mix powder and liquid manually in a clean Multivac (evacuation unit) mixing bowl using a spatula until no more dry powder is present in the bowl (no more than 30 seconds). Alternatively, use the Multivac compact mixing unit for pre-mixing. Close the mixing bowl. Place the bowl in an evacuation unit until the maximum vacuum is reached. Then stir the investment compound under vacuum for 60 s.

For large-span bridges (6 or more units) or complicated constructions, the mixing time in vacuum can be shortened to 30 seconds if required.

**CAUTION!** If you use the Multivac compact for mixing, use Program #1, as for **Deguvest SR**, or the standard program. Do not use the Optivest program.

### 4.3 Pouring silicone duplicating moulds

Depending on the cast, 100 –150 g powder are needed (**100 g : 15 ml or 150 g : 22.5 ml**). It is easier to pour silicone duplicating casts if a wetting agent is first applied. **Caution:** The moulds must be blow-dried completely before pouring. Silicone duplicates should **not** be poured immediately after duplicating the cast. Removing the cast creates a deformation that takes some time to disappear. Allow the duplicate to stand for at least 15 minutes before pouring the **Optivest SR**.

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the duplicating mould has been filled, turn off the vibrator immediately. Allow to harden for 30 minutes. Carefully remove the cast from the duplicating mould.

### 4.4 Pouring agar-agar duplicating moulds

Water the cast in hand-warm water for approximately 10 minutes before applying the duplicating gel. After watering, blot the liquid carefully off the cast or blow-dry the cast. For gel duplication, follow the gel manufacturer's instructions.

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the duplicating mould has been filled, turn off the vibrator immediately. Allow to harden for 30 minutes. Carefully remove the cast from the duplicating mould!

#### 4.5 Preparation and drying

To create undercuts that improve the connection between the cast and the jacket, carefully bevel the bottom (base) edges of the refractory cast using a dry dental trimmer.

Dry the cast in a drying chamber or oven pre-heated to approximately 200 °C for about 40 minutes. Allow the cast to slowly cool to room temperature on the bench. Correctly dried casts exhibit a light colour and absorb bonding agents well.

#### 4.6 Hardening/Preparation

When using silicone duplicating moulds, **Optivest SR** refractory casts do not have to be hardened. A bonding agent can be used to achieve better adhesion of the wax model on the cast. Make sure the bonding agent does not form a layer of its own on the surface of the cast, to avoid the formation of casting fins, filling up of retentive recession and compromising fit. Check the consistency of the bonding agent from time to time.

If the refractory cast was produced using duplicating gel, hardening by immersion in a hardener bath is recommended. Make sure to follow the pertinent instructions for these products. Old hardener can be viscose and touch and form an unwanted layer on the model.

**No** expansion liquid must be brushed onto the cast (or applied by immersion) on the surface of the cast after drying. Do not start modellation until the model has cooled down to hand temperature!

## 5 Creating the casting mould

### 5.1 Stress relief

In most cases, the wax model does not require stress relief. If the wax model has very delicate details, wax objects may be painted or sprayed with a thin layer of wetting agent and then carefully blown dry (do not use pressurized air). The wax surface must be absolutely dry before investing.

### 5.2 Protective investment coating

Before investing the wax model in **Optivest SR**, the wax surface can be pre-treated with a protective investment coating to ensure that no bubbles are formed in undercuts or in complicated areas of the casting, replicating surfaces exactly as they were on the wax model.

### 5.3 Types of flask

We recommend using reusable plastic or silicone flasks. These ensure an even distribution of the investment compound around the model and a uniform distance between the base of the cast and the outer surface of the investment. This distance should be at least 10 mm. It is important for the stability of the muffle during the heating phase and at the same time reduces the risk of muffle tears.

### 5.4 Mixing the investment compound

Measure **Optivest SR** and mixing liquid and mix as described in 4.2.

For mixing ratios refer to the recommended concentrations in the table following 6, Expansion control.

### 5.5 Investing

Pour mixed investment compound into the muffle using slight vibration. Make sure to void trapping air. After the casting mould has been filled, turn off the vibrator immediately.

When investing crown and bridge frameworks, it is recommended to fill the crowns with investment first. Optivest SR is primarily an investment intended for metal casting. By contrast with specialized crown and bridge investment compounds, the working time is relatively short.

### 5.6 Hardening

Allow the casting mould to harden for approximately **45 minutes**. Once the investment compound has set, slightly grind the cover of the muffle using a dry trimmer to remove the glaze layer so that any gas can escape.

### 5.7 Wax elimination

Pre-heat the furnace to 250 °C and place the muffle into the furnace with its funnel facing down. Keep the muffle in the chamber for 45–60 minutes to eliminate the wax! Make sure to keep the chamber door closed to keep noxious gases from escaping!

### 5.8 Pre-heating

Once the wax has been eliminated, heat the mould to the final temperature for metal casting (1000–1050 °C). The heating rate should not exceed 8 °C/min. The final temperature should be held for 60 minutes (longer if a larger number of muffles is used – add 10 minutes per muffle). Muffles can also be pre-heated overnight, under program control.

### 5.9 Casting

Follow the instructions of the alloy manufacturer and the instructions for your casting device.

### 5.10 Divesting

Allow the flask to cool to room temperature after casting (for at least 30 minutes).

**Caution:** After pre-heating, **Optivest SR** will lose its low-dust properties! We recommend rinsing the muffle approximately 15 minutes before divesting to reduce the amount of dust released.

Never quench the muffle in cold water. Doing so will damage the metal alloy. Cleave it repeatedly with a pair of plaster pliers or a divesting chisel and remove from casting. The framework can then be sandblasted with a suitable abrasive (150-250 µm).

## 6 Expansion control

Table 1: Mixing

Mixing ratio (powder : liquid)	100 g : 15 ml
Mixing time (under vacuum)	30 – 60 s
Multivac compact	Program #1 (Deguvest SR)

We have summarized the most important information for you in a table. We would like to emphasize that these are only **recommendations** developed under controlled laboratory conditions. Variations in temperature ranges, modeling materials, procedures, muffle systems, pre-heating ovens, casting units, and many other factors may affect the results. If the results you achieve are unsatisfactory, consult your technical adviser, who will be happy to assist you.

Table 2: Metal casting

Form	Material	Powder amount (g)	Concentration (%)	Standard CoCrMo alloys (Biosil)		Concentration (%)	High-gold casting alloys (Degunorm classic, Degulor M)	
				ml liquid	ml distilled water		ml liquid	ml distilled water
Refractory cast	Silicone	100	75	11	4	75	11	4
		200	75	22	8	75	22	8
	Gel	100	65	9	6	65	9	6
		200	65	18	12	65	18	12
Investment cover	Resin muffle former	400	Like the cast			Like the cast		
		450	Like the cast			Like the cast		
	Crepe sleeve	400	Like the cast			Like the cast		
		450	Like the cast			Like the cast		

Table 3: Single and multi-unit restorations

		Muffle size	Concentration (%)	Pd-based alloys (DegupalG, DegustarF)		Concentration (%)	AgPd alloys (Pallig M)	
				ml liquid	ml distilled water		ml liquid	ml distilled water
Inlays, Onlays, partial crowns	Wax	1x	50	7	7	50		
		3x	40	9	13	40		
		6x	40	18	27	40		
Crowns and bridges Primary crowns	Wax	1x	50	7	7	50	7	7
		3x	40	9	13	40	9	13
		6x	40	18	27	40	18	27
		9x	30	20	47	30	20	47
	Swaging foil	1x	40	6	9	40	5	9
		3x	30	6	16	30	7	15
		6x	30	13	32	30	14	31
Secondary crowns	Wax	1x	60	9	6			
		3x	60	18	12			
		6x	60	27	18			
		9x	50	34	33			
	Modelling resin	1x	80	12	3	80		
		3x	70	16	6	70		
		6x	70	31	14	70		
		9x	60	40	27	60		

Table 4: Non-precious alloys

**Optivest SR is compatible with all commercial non-precious alloys (with the exception of titanium). Follow the alloy manufacturers' recommendations for pre-heating temperatures. The concentrations listed here are recommendations only. Because of the many different alloys and their widely divergent properties, we recommend reference castings in the case of complicated objects.**

		Non-precious alloys	
Muffle size	Concentration (%)	ml liquid	ml distilled water
1x	80	12	3
3x	80	19	4
6x	90	41	4
9x	90	60	7
1x	80	12	3
3x	80	19	4
6x	90	41	4
9x	90	60	7

We wish you much success with **Optivest SR** investment. If you have any questions as regards handling or technical details, please do not hesitate to call your local DeguDent representative, or visit our website at [www.degudent.com](http://www.degudent.com).

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