Phantom-Metall

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Instructions for use

Mode d’emploi

Istruzioni d’uso

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DeguDent GmbH

A. Beschreibung

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A. Description

Phantom-metal is a dental training material which is in its processing and handling very close to the characteristics of commercial dental precious-metal alloys. **Caution: phantom-metal is not mouth resistant!** For this reason it should be used only to carry out practice exercises and to produce illustrative or demonstration models.

Phantom-metal is supplied in the form of ingots. In addition, two types of solder matched to the alloys are available:

Phantom-metal solder 1 (main solder)
Phantom-metal solder 2 (secondary solder)

B. Technical Data

Melting range °C	Hardness HV 5	Tensile strength * N/mm²	Yield strength * N/mm²	Elongation at fracture* %	Density g/cm³
900–1030	130	400	230	50	9.0

*Test specimens in accordance with DIN 13 906

Phantom-metal solder 1	Working temperature 710 °C
Phantom-metal solder 2	Working temperature 660 °C

C. Instructions for Use

- Embed the wax model in a commercial investment such as, for example, Deguvest® HFG, Deguvest® F or Deguvest® CF.
- Drying time of the casting ring, depending on the size of the ring, between 30 and 45 min.
- Time required for wax elimination, depending on the size of the casting ring, between 30 and 60 min at 300 °C.
- Preheat the casting rings to 700 °C. Maintain this final temperature for 20 to 60 min depending on the size of the ring.
- Melting
 - a) with the electric high heat melting device associated with the e.g. Multicas® compact or Prestomat® compact:
 - Preheat the graphite crucible and fill it with ingots. In order to prevent the formation of slag, the addition of Veriflux® melting powder is advisable.
 - For melting purposes, set the control unit in the melting system to 1050–1100 °C. After the alloy has been molten, superheat for 30 sec;
 - b) with the propane/oxygen flame in the universal crucible:
 - Sprinkle Veriflux® melting powder on the ingots and melt with a reducing flame. After the alloy has been molten, superheat for 5 sec. longer.
- Soldering
 - Grind the surfaces to be soldered so that they are metallically bright. Assemble soldering parts that are as narrow as possible. In order to solder several bridge parts, after they have been fixed in position embed them in a soldering block made of a soldering investment such as, for example, Deguvest® L. After the adhesive wax has been boiled off, coat the soldering gaps and their vicinity, as well as the phantom-metal solder, with Anoxan® flux. Solder as usual. Remove the flux residues.
- Pickling
 - Phantom-metal can be pickled **only conditionally** in NeaciD® or in dilute sulfuric acid. Therefore we recommend that the surface of the cast objects be cleaned mechanically by lapping with emery powder, rubberizing, and polishing.

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