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Anoxan[®]

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Xn

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Anoxan®

Anoxan is a pasteous fluxing agent used for soldering in crown and bridge procedures.

It is suitable for solders with a working temperature up to approx. 850 °C.

Note: Anoxan is not suitable for soldering procedures following ceramic firing.

Anoxan is a water-based product, making the burning off of a fatty base along with the associated formation of flames unnecessary.

Anoxan dissolves existing layers of oxides and prevents formation of new ones.

Anoxan causes excellent flow of solders and does not bubble during heating.

Anoxan is economical in use as it contains the main constituents in a concentrated form.

Availability: Jars of 50 ml. each (approx 5 oz).

REF 53 2529 0001

Instructions for use:

1. Clean solder joints and adjacent areas thoroughly.
2. Join the parts to be soldered without inducing any stress by positioning with wax or with pattern resin and invest using Deguvest® L soldering investment material.
3. Boil out wax or burn out the resin in the preheating furnace and apply Anoxan to exposed parts to be soldered as well as the solder.
4. Moisten also the solder before positioning on the objects.
5. Solder as usual.
6. Remove remains of Anoxan in hot Neacid® solution.
7. Close Anoxan container tightly after use and prevent from drying.

Attention!

Harmful if swallowed. Keep out of reach of children. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable gloves and eye/face protection.

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DeguDent GmbH