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# Biosint<sup>®</sup> Extra

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso

**DeguDent**

*A Dentsply International Company*

# Biosint® Extra

<b>D</b>	Gebrauchsanweisung .....	3
<b>GB</b>	Instructions for use .....	8
<b>F</b>	Mode d'emploi .....	13
<b>I</b>	Istruzioni per l'uso .....	18
<b>E</b>	Instrucciones para el uso .....	23

# Biosint® Extra

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# Biosint® Extra

Biosint® Extra is a gypsum-free, phosphate-bonded special investment compound for the production of duplicate casts and dies for the CoCr model cast technique.

Biosint Extra is mixed with Biosol® E.

Adding the special Biosol E mixing liquid to water changes the total expansion.

Mixing only with water (distilled water) is possible.

## 1 Safety information

Mixing with other liquids or materials during the use must not be carried out, since during this process health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

Carry out the “coarse” removal of the casting after the casting operation under water or soak the investment material mold.

## 1.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

## 2 Technical data

Mixing ratio	100 g powder : 15 ml liquid
Total expansion	0.95 % – 1.65 %
Mixing time under vacuum	60 sec
Processing time span	3 – 5 min
Compression strength	according to concentration of the expansion liquid between 15 and 20 MPa

The total expansion, i.e. the sum of setting expansion and thermal expansion can be changed by diluting the Biosol E mixing liquid.

Mixing liquid consisting of:		Concentration	Total sum of:	
Water	Biosol E	Biosol E	setting exp.	therm. exp.
4 parts	0 part	0 %	0.30 %	0.65 %
2 parts	2 parts	50 %	0.35 %	0.95 %
1 part	3 parts	75 %	0.45 %	1.05 %
0 part	4 parts	100 %	0.50 %	1.15 %

### 3 Important information

The safety information mentioned under 1 and 1.1 must be observed prior to the usage.

Biosint Extra must not be brought into contact with gypsum or gypsum-containing materials since this may lead to sulphur damage of high-melting alloys. Therefore special mixing cups and spatulas must be used for mixing Biosint Extra.

Biosint Extra must be stored under cool and dry conditions.

If the product is stored properly, the storage capacity of the investment material amounts to 18 months.

Biosil E must be stored at a room temperature between 5 °C and max. 30 °C.

### 4 Production of duplicating models

#### 4.1 Mixing of the investment compound

Premix the powder and liquid by hand with a spatula in the clean, freshly rinsed cup of the Multivac® evacuating unit. Then let the investment compound stand under vacuum for a short time. Subsequently mix under vacuum for approx. 60 sec.

#### 4.2 Pouring the duplicating mould

Mix the investment compound flow into the duplicating mould under slight vibration and carefully avoid the formation of bubbles. After the duplicating mould is filled, immediately switch off the vibrator.

**Note:** Silicone duplicating mould can be poured easier if they have been wetted previously with Waxit®. The mould must be blown dry before pouring with investment compound.

Investment compound models which have to be produced in silicone duplicating moulds must be mixed with higher concentration of Biosol E; they require greater expansion.

### 4.3 Setting

Let the investment compound set for approx. 30 min and then remove carefully from the duplicating mould.

### 4.4 Drying

Dry the model for approx. 40 min in a drying furnace preheated to around 220 °C.

### 4.5 Hardening

**Note:** Investment compound models which have been produced in silicone duplicating moulds need not to be hardened.

Immerse the dried investment compound model in cold immersion hardener for approx. 10 s.

Let the immersion liquid drip off the cast and put the cast back in the switched-off drying furnace for approx. 5 min.

Using hardening waxes, hardening sprays or wax adhesives is also possible.

### 4.6 Modelling

Modelling on hand-warm model.

## 5 Production of the casting mould

### 5.1 Relaxing

Thinly brush the wax object with the Waxit wetting agent and then carefully blow it dry.

### 5.2 Fine investment

Before investing of the wax modelation with Biosint Extra, the entire wax surface can be pretreated with fine investment compound. This provides the best guarantee of avoiding air bubbles at undercuts or points on the poured object which are difficult to see and also results in exact surface reproduction of the cast.

### 5.3 Mixing of the investment compound

Dose Biosint Extra and mixing liquid and mix as described under 4.1.

We recommend as mixing ratio 100 g powder : 15 ml liquid.

**Note:** The casting mould should be produced in the same or somewhat lower concentration as the model!

### 5.4 Investing

Let the mixed investment compound flow into the mould form under light vibration and carefully avoid the formation of bubbles.

Immediately switch off the vibrator after the mould is filled.

### 5.5 Setting

Let the mould set for approx. 30 min.

### 5.6 Burn out of the wax

Place the mould in a furnace preheated to 300 °C and burn out the wax for approx. 45 – 60 min.

### 5.7 Preheating

Then bring the mould to the required final temperature of 1000 – 1050 °C required for CoCrMo alloys. Maintain the final temperature for approx. 60 min.

The notes provided under 5.6 and 5.7 also apply accordingly for program-controlled preheating overnight.

### 5.8 Casting

Cast according to the processing instructions of the alloy manufacturer.

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**DeguDent GmbH**