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


# Biosint<sup>®</sup> Supra

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso  
Инструкция по  
применению



**DeguDent**

*A Dentsply International Company*

# Biosint®-Supra

<b>D</b>	Gebrauchsanweisung .....	5
<b>GB</b>	Instructions for use .....	10
<b>F</b>	Mode d'emploi .....	15
<b>I</b>	Istruzioni per l'uso .....	20
<b>E</b>	Instrucciones para el uso .....	25
<b>RUS</b>	Инструкция по применению .....	30

D

# Biosint®-Supra

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# Biosint®-Supra

Biosint®-Supra is a plasterfree, phosphate-bonded special investment compound for the production of duplicate casts and dies for the CoCr model cast technique.

Biosint-Supra is mixed with Biosol® E.

Adding water (aqua dest.) to Biosol E mixing liquid changes the total expansion.

Mixing only with water (distilled water) is possible.

Biosol E must be ordered separately.

## 1 Technical data

Mixing ratio	100 g powder : 15 ml liquid
Total expansion	0.95 – 1.65 %
Mixing time under vacuum	60 sec
Working time range	3 – 5 min
Compressive strength	between 15 and 20 MPa (depending on the concentration of the mixing liquid)

The total expansion, i.e. the sum of setting expansion and thermal expansion can be changed by diluting the Biosol E mixing liquid.

Mixing liquid consisting of:		Concentration	Total sum of:	
Water	Biosol E	Biosol E	setting exp.	therm. exp.
4 parts	0 part	0 %	0.30 %	0.65 %
2 parts	2 parts	50 %	0.35 %	0.95 %
1 part	3 parts	75 %	0.45 %	1.05 %
0 part	4 parts	100 %	0.50 %	1.15 %

## 2 Important information

The safety information mentioned under 5. and 5.1 must be observed prior to usage.

Biosint-Supra should not be brought into contact with plasters or plaster-containing materials since this may lead to sulphur damage of high-melting alloys. Therefore separate mixing bowls and spatulas must be used for mixing Biosint-Supra.

Biosint-Supra must be stored under cool and dry conditions.

If the product is stored properly, the storage capacity of the investment material amounts to 18 months.

Biosil E must be stored at a room temperature between 5 °C and max. 30 °C.

## 3 Production of duplicating models

### 3.1 Mixing of the investment compound

Premix the powder and liquid by hand with a spatula in a clean, freshly rinsed bowl of the Multivac evacuating unit. Then let the investment compound stand under vacuum for a short time. Subsequently mix under vacuum for approx. 60 sec.

### 3.2 Pouring the duplicating mould

Refractory models made using Agar-Agar gel should be mixed with 50 % liquid concentration.

Let the mixed investment compound flow into the duplicating mould under slight vibration and carefully avoid the formation of bubbles. After the duplicating mould is filled, immediately switch off the vibrator.

**Note:** Silicone duplicating moulds can be poured easier if they have been wetted previously with Waxit®. The mould must be blown dry before pouring with investment compound.

Investment compound models which are made in silicone duplicating moulds must be mixed with a higher concentration of Biosol E – we recommend a concentration of 60 %.

### 3.3 Setting

Let the investment compound set for approx. 30 min and then remove carefully from the duplicating mould.

### 3.4 Drying

Dry the model for approx. 40 min in a drying furnace preheated to around 220 °C.

### 3.5 Hardening

**Note:** Investment compound models made in silicone moulds need not be hardened.

Immerse the dried investment compound model in the Sintadur cold immersion hardener for approx. 10 s.

Let the immersion liquid drip off the cast and put the cast back in the switched-off drying furnace for approx. 5 min.

Use of hardening waxes, hardening sprays or wax adhesives is also possible.

### 3.6 Modelling

Begin modelling on hand-warm model.

## 4 Production of the casting mould

### 4.1 De-greasing

Thinly brush the wax object with the Waxit wetting agent and carefully blow it dry.

### 4.2 Fine investment

Before investing the wax modelation with Biosint-Supra, the entire wax surface can be pretreated with fine investment compound. This gives the best guarantee of avoiding air bubbles at undercuts or points on the poured object which are difficult to see and also results in exact surface reproduction of the cast.

#### 4.3 Mixing the investment compound

Dose Biosint-Supra and mixing liquid and mix as described under 3.1.

We recommend as mixing ratio 100 g powder : 15 ml liquid.

**Note:** The casting mould should be made with the same or somewhat lower concentration as the model!

#### 4.4 Investing

Let the mixed investment compound flow into the mould form under light vibration and carefully avoid the formation of bubbles.

Immediately switch off the vibrator when the mould is full.

#### 4.5 Setting

Let the mould set for approx. 30 min.

#### 4.6 Burn out of the wax

Place the mould in a furnace preheated to 300 °C and burn out the wax for approx. 45 – 60 min.

#### 4.7 Preheating

Heat the mould to the final temperature of 1000 – 1050 °C required for CoCrMo alloys. Maintain the final temperature for approx. 60 min.

The notes given under 4.6 and 4.7 also apply accordingly for program-controlled preheating overnight.

#### 4.8 Casting

Cast according to the processing instructions of the alloy manufacturer.

## 5 Safety information

Mixing with other liquids or materials during this process means health hazards, as well as impairments to the casting quality, cannot be excluded.

All dental investment compounds contain crystalline quartz and/or cristobalite.

These substances can adversely affect health if inhaled. Can cause silicosis.  
Always use a sufficiently powerful extractor or wear a mask.  
Carry out the "coarse" removal of the casting under water or soak the investment material mold beforehand.

### **5.1 Information for the processing of mixing liquids for dental investment materials**

Following contact with skin:

Wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

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