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


Cergo fit SPEED

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Instruções de utilização
**Инструкция по
применению**




DeguDent

A Dentsply International Company

Cergo fit SPEED

D	Gebrauchsanweisung	5
GB	Instructions for use	8
F	Mode d'emploi	11
I	Istruzioni per l'uso	14
E	Instrucciones para el uso	17
P	Instruções de utilização	20
RUS	Инструкция по применению	23

Cergo fit SPEED

D

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Cergo fit SPEED

Cergo fit SPEED (REF 2532 0200) is a phosphate-bound, rapid burnout investment, for pressable ceramics. Cergo fit SPEED has no plaster or graphite. Cergo fit SPEED is mixed with Cergo fit SPEED liquid (REF 2540 0117).

1 Processing Recommendations

Processing Parameters	Recommended
Material Temperature	21 °C / 70 °F
Mixing Ratio	100 g : 25 ml
Vacuum Mixing Time	1 min
Working Time	5 – 7 min

2 Important Tips

Please strictly observe the safety tips in Paragraph 8 before use. Cergo fit SPEED must not come into contact with plasters or compounds containing plaster, because this can cause damage to the press ceramics (Cergo). Separate mixing bowls and clean spatulas must therefore be used for mixing Cergo fit SPEED.

Store Cergo fit SPEED in a cool and dry place. If stored properly, investment will keep for 18 months.

Fully protect Cergo fit SPEED liquid from frost, and store at a room temperature between 5 °C/41 °F, and no more than 30 °C/86 °F!

3 Expansion Control

Cergo fit SPEED is mixed with Cergo fit SPEED liquid. Expansion can be controlled by diluting this liquid with distilled water, thus achieving an optimal fit for any geometry.

Specific concentration *recommendations* to be individually adapted:

- Incisor restoration: approx. 80 % → 20: 5 (ml)
- Premolar restoration: approx. 75 % → 19: 6 (ml)
- Inlay, triple surface: approx. 85 % → 21: 4 (ml)
- Inlay, double surface: approx. 90 % → 22: 3 (ml)
- Inlay, single surface: approx. 90 % → 22: 3 (ml)
- Veneer: approx. 85 % → 21: 4 (ml)

Please note that many processing parameters (e.g. temperature, modelling material etc.) can influence the fit and deviations from the aforementioned recommendations are therefore possible.

4 Mixing the Cergo fit SPEED Investment

Premix powder and liquid in a clean mixing bowl of the Multivac® stirring and evacuation unit for about 30 seconds (including period of sprinkling by hand) with a clean spatula. Let the investment stand under vacuum for about 30 sec. Next, mix under vacuum for about 60 sec. with the Multivac.

When mixing with the Multivac compact, please use Program No. 1 (Deguvest CF).

When mixing Cergo fit SPEED Investment, please observe the following:

The mixing ratio and mixing time are important for consistent results. We recommend that the mixing ratios and times provided be complied with exactly. High vacuum in the mixing bowl determines quality; only start the mixing procedure after waiting for a full vacuum display.

4.1 Investing

Place ring former on sprue base former. While gently vibrating ring assembly, pour the vacuum-mixed investment into the ring, allowing the investment to flow bubble free into and over wax patterns. Place the ring leveller device into the ring former. The excess investment will be squished up and through the opening in the leveller device.

Placement of the leveller insures correct ring height, and proper alignment of the plunger and the ring.

After **10 minutes**, the ring may be removed from the ring former. Be sure to remove the uneven nub from end of ring with knife, making sure you have a flat and even surface!

5 Burnout

Cergo fit SPEED is suitable for so-called rapid burnout techniques. In order to achieve consistent results in this heating procedure, please observe the following:

After **15 min.**, measured from the moment the powder and liquid are mixed together, the ring must be placed immediately into the furnace preheated to **850 °C**. Lean the ring against the back or side walls of the oven at a **45°** angle.

In order to achieve homogenous heating of the ring, the holding time of **60 min.** must be adhered to at the final temperature.

6 Pressing

Pressing is done in accordance with the processing instructions of the particular press ceramic manufacturer.

7 Divesting

After pressing, let the ring cool in the open air at room temperature. Stand the cooled ring on its end on the plunger. Place a plunger next to it and lay the tip of a pencil over the edge slightly. Rotate the ring, engaging the pencil lead to the side of the ring, scoring a line all around. Use a cut-off disk at the pencil line, and cut off plunger from ring. Also, cut the investment off of the plunger.

Sandblast the ring with **4 bar** pressure until the units become visible, then drop the pressure to about **2 bars**. Use only 50 micron glass bead grit (**50 µm**).

8 Safety Tips

Do not mix with other liquids or materials during use, because this can lead to health risks or reduced compression quality.

All dental investments contain portions of crystalline quartz and/or cristobalite. These substances may cause damage to health and/or silicosis. In general, only work with respiratory protection.

8.1 Information for the processing of mixing liquids for dental investment compounds

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

Further information is contained in the safety data sheet.

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DeguDent GmbH