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Deguvest[®] GF

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Deguvest® GF

| | | |
|-----------|-------------------------------------|----|
| D | Gebrauchsanweisung | 5 |
| GB | Instructions for use | 11 |
| F | Mode d'emploi | 17 |
| I | Istruzioni per l'uso | 23 |
| E | Instrucciones para el uso | 29 |

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Deguvest GF is a carbon-free, phosphate-bonded and plasterfree, special investment especially well suited for the precision casting of carbon-sensitive, precious metal alloys.

1 Technical data

Deguvest GF

| | |
|--|---------------|
| Mixing ratio powder : liquid | 100 g : 22 ml |
| Total linear expansion | 1.2 – 2.4 % |
| Mixing time under vacuum | 60 sec |
| Working time range | 4 – 6 min |
| Compressive strength depending on the concentration of the mixing liquid | 4 – 8 MPa |

| Mixing liquid composed of | | Concentration | Setting expansion | Thermal expansion |
|---------------------------|------------------------|---------------|-------------------|-------------------|
| Distilled water | Deguvest mixing liquid | | | |
| 3 parts | 1 part | 25 % | 0.30 % | 0.90 % |
| 2 parts | 2 parts | 50 % | 0.60 % | 1.00 % |
| 1 part | 3 parts | 75 % | 0.90 % | 1.10 % |
| | 4 parts | 100 % | 1.15 % | 1.25 % |

Setting expansion + thermal expansion (mean values) = total expansion

2 Important information

The safety information mentioned under 6 and 6.1 must be observed prior to use.

Deguvest GF must not get in contact with plasters or plaster-containing materials since this can cause sulphur damage of high-melting alloys. Therefore special mixing bowls and spatulas have to be used for mixing Deguvest.

Deguvest must be stored under cool and dry conditions.

If the product is stored properly, the storage time of the investment material amounts 18 months.

Deguvest mixing liquid must be stored at a room temperature between 5 °C and max. 30 °C.

Typical concentrations of the mixing liquid for investing of objects in Deguvest GF

| Object: | MOD-Inlay | anterior crown | anterior crown | posterior crown | posterior crown | 6° tapered secondary | 6° tapered secondary | telescopic secondary | telescopic secondary | RS-circular | bridge with 8 members |
|---|-----------|----------------|----------------|-----------------|-----------------|----------------------|----------------------|----------------------|----------------------|-------------|-----------------------|
| Modellation material: | wax | wax | Adapta | wax | Adapta | Adapta | pattern | Adapta | pattern | pattern | wax |
| Alloys: examples in () | | | | | | | | | | | |
| High gold content c & b alloys (Degulor M) | 60–70% | 65–75% | 65–75% | 65–75% | 65–75% | 50–60% | 40–50% | 60–70% | 55–65% | 50–60% | 65–75% |
| Reduced c & b alloys (Stabilor NF IV) | 60–70% | 65–75% | 65–75% | 65–75% | 65–75% | 50–60% | 40–50% | 60–70% | 55–65% | 50–60% | 65–75% |
| High gold content ceramic bonding alloys (Degudent U) | 60–70% | 65–75% | 65–75% | 65–75% | 65–75% | 50–60% | 40–50% | 60–70% | 55–65% | 50–60% | 65–75% |
| Reduced ceramic bonding alloys (Deva 4) | – | 65–75% | 65–75% | 65–75% | 65–75% | – | – | – | – | 55–65% | 70–80% |
| Palladium-based alloys (Degupal G) | – | 70–80% | 70–80% | 70–80% | 70–80% | – | – | – | – | 60–70% | 75–85% |

3 Expansion control

Deguvest investments are mixed with Deguvest mixing liquid (does not apply to Deguvest California). To control the expansion of Deguvest the liquid has to be diluted with distilled water. By this the concentration of the respective alloy can be compensated.

For detailed application recommendation see table on page 12.

4 Manufacturing of the casting mold

- Dilute Deguvest mixing liquid to the desired concentration in a mixing bottle by adding **distilled** water.
- Line the casting ring using a moistened casting ring liner or Deguvest fleece.
- Fix the wax modellation in place with wax on the sprue former.

- If you use a wetting agent then ensure the surfaces of your copings are free from any residue before investing.
- Measure the Deguvest powder and mixing liquid in the following ratios:

100 g powder to 22 ml mixing liquid

or corresponding to the data in the following table.

| Casting ring | Deguvest GF -powder | Deguvest- mixing liquid |
|--------------|------------------------|----------------------------|
| 1 x | 1 x 60 g = 60 g | 13 ml |
| 3 x | 1 x 150 g = 150 g | 33 ml |
| 6 x | 2 x 150 g = 300 g | 68 ml |
| 9 x | 3 x 150 g = 450 g | 102 ml |

4.1 Mixing the investing compound

Premix powder and liquid for approx. 30 sec by hand with a spatula in a clean bowl of the Multivac® stirring and evacuating unit. Let the investing compound stand for about 30 sec under vacuum. Finally, mix for about 60 sec under vacuum with the Multivac®.

During mixing of the investing compound, the following points, among others, must be observed:

- High degree of vacuum in the mixing bowl determines quality; wait until the full vacuum is indicated, and only then start the mixing process.
- The mixing ratio influences the material characteristics.
A mixing ratio outside the recommendation has the following effects:
a larger amount of liquid added has the following effects:
 - thinner consistency
 - longer working time range
 - slightly lower setting expansion
 - slightly lower strength
 a smaller amount of liquid added has the following effects:
 - smoother surfaces of the cast
 - somewhat higher expansion

- The mixing time is important for consistent results; observe the recommended mixing time exactly.

Mixing time shorter:

- uncontrollable setting expansion and strength
- air bubbles can remain in the investing compound and result in casting beads

Mixing time longer:

- lower setting expansion
- shorter setting times

4.2 Investing

Let the mixed investment material flow into the casting ring under light vibration and carefully avoid the formation of bubbles. Immediately turn off the vibrator after the casting ring has been filled.

4.3 Setting time

Let the casting mold set for 20 – 50 min (see table, Heat treatment).

5 Preheating and Casting

5.1 Wax elimination

Place the casting mold in a furnace preheated to about 300 °C, and eliminate the wax in about 45 – 60 min (see table, Heat treatment).

In order not to impair the result by wax residues in the casting mold, the wax should be eliminated immediately after the investment material has set.

5.2 Preheating

Heat the casting mold to the final temperature required for the particular alloy. Hold the final temperature for about 20 – 60 min (see table, Heat treatment).

A factor that influences the casting result is the uniform heating of the casting ring to the alloy-specific preheating temperature. In order to make optimum use of the thermal expansion and to avoid the formation of cracks, the heating should be carried out slowly, especially in the ranges of the cristobalite and quartz phase shifts.

The instructions given above also apply analogously for the program-controlled preheating.

5.3 Preheating with holding times

The outstanding accuracy of fit and surface quality of the Deguvest investment compounds can not be improved by preheating with several holding times (additional “stops” before reaching the final temperature).

A uniform increase of temperature is important; in the ranges of the cristobalite and quartz phase shifts (at about 270 °C and 560 °C respectively) the increase in temperature should be carried out slowly (ideally 7 °C/min).

The parameters that influence the casting result, such as modelling material (wax, resin, adapter cap), geometrical form and size of the dental object, casting ring size, and alloy must be regarded individually with each case, and will have a great influence on the quality of the casting.

– Holding times do not change these factors of influence!

| Heat treatment | Size of the casting ring | | | |
|--|--------------------------|--------|--------|--------|
| | 1 x | 3 x | 6 x | 9 x |
| Setting/drying time | 20 min | 30 min | 40 min | 50 min |
| Wax elimination temperature | 300 °C ** | | | |
| Duration of wax elimination | 30 min | 40 min | 50 min | 60 min |
| Initial temperature of the electric preheating furnace | 300 °C ** | | | |
| Preheating time after final temperature is reached * | 20 min | 30 min | 45 min | 60 min |

* alloy-specific

** When the casting rings are placed in a preheating furnace preheated to 300 °C, a slow distribution of heat takes place inside the casting mold because of the thermal conductivity of the casting mold

5.4 Casting

Cast according to the directions for use from the manufacturer of the alloy.

5.5 Devesting

After casting let the casting ring cool down to room temperature and devest. For this purpose the investment material mold must be soaked in water for approx. 15 min and carefully split up at several points using plaster pliers. So the dental casting can be easily removed from the investing compound, and formation of dust is prevented. Final residues of investing compound in the crowns can be picked out or carefully sandblasted with aluminum oxide, grain size 100 to 150 µm.

Do not use a hammer for devesting!

6 Safety information

Mixing with other liquids or materials during this process health hazards as well as impairments to the casting quality cannot be excluded.

All dental investment compounds contain parts of crystalline quartz and/or cristobalite.

These substances can adversely affect health if they are breathed in. Always use a sufficiently powerful extractor or wear a mask.

Carry out the “coarse” removal of the casting after the casting operation under water or soak the investment material mold.

6.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

For further details, please refer to the Safety Data Sheet.

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