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Deguvest[®] Impact

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DeguDent
A Dentsply Company

Deguvest[®] Impact

D	Gebrauchsanweisung	5
GB	Instructions for use	11
F	Mode d'emploi	17
I	Istruzioni per l'uso	23
E	Instrucciones para el uso	29

Deguvest[®] Impact

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Deguvest[®] Impact

Deguvest Impact is a carbon-free, phosphate-bonded and plasterless universal investment material which is suitable for the entire spectrum of precious metal precision casting.

Deguvest may be used for the conventional as well as for the rapid burnout method heating. For this purpose the special liquid must be used. Carefully read the information provided under 6.3.

1 Safety information

Mixing with other liquids or materials during processing health hazards as well as impairments to the casting quality cannot be excluded.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

1.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

2 Technical data

	Deguvest Impact
Mixing ratio powder : liquid	100 g : 23 ml
Total linear expansion	1.2 – 2.4 %
Mixing time under vacuum	60 s
Processing time	4 – 6 min
Compressive strength (depending on concentration of the mixing liquid)	4 – 8 MPa

Mixing liquid composed of		Concentration	Setting expansion	Thermal expansion
Distilled water	Deguvest mixing liquid			
3 parts	1 part	25 %	0.30 %	0.90 %
2 parts	2 parts	50 %	0.60 %	1.00 %
1 part	3 parts	75 %	0.90 %	1.10 %
	4 parts	100 %	1.15 %	1.25 %

Setting expansion + thermal expansion (mean values) = total expansion

3 Important information

Be sure to observe information concerning safety prior to use under sections 1 and 1.1.

Deguvest Impact must not come into contact with plasters or plaster-containing materials since this may lead to sulphur damage of heat melting-alloys. Therefore special mixing bowls and spatulas must be used for mixing Deguvest Impact.

3.1 Investment Storage

Most of you are probably aware of the problems which can occur when investment is subjected to variations in environment temperature (e.g. summer and winter) and the effect this can have on your castings.

That is why it is very important to find the right place to store your supply of investment. The investment should be stored in a cool and dry place, between 20 °C – 23 °C (71.6 °F – 73.4 °F). If you have a place in your laboratory where these conditions can be found throughout the year, then that's the place to store your Investment!

If it is not possible for you to store your investment at this temperature, then you should at least make sure that the work bench where investing is done fits these requirements (use a temperature control cupboard or similar).

Please note the symbols on product labels:

REF Product code

LOT Batch or lot number

 Use before

 Consult instructions for use

Avoid any nasty surprises by observing the following:

- Avoid using any investment that has just been brought out of an especially warm or cold place. (e.g. storage rooms, silos, cupboards etc.)
- Fill up stocks at your workbench before your supply runs out.
- Give the investment time to adjust to the environment!
- Avoid using investment that has just been delivered to the laboratory. Again, give the investment time to acquire the right temperature.

- Sporadic temperature checks at the workbench help maintain constant temperatures.
- If space allows, stock up on investment (especially LIQUID) before the cold winter months arrive.
- This way, you forego the risk of exposing the liquid to frost during delivery.

Typical concentrations of Deguvest® Liquid for investing objects in Deguvest® Impact

Please note that the special liquid must be used for rapid heating (see 6.3).

Object:	MOD-Inlay	Anterior crown	Anterior crown	Posterior crown	Posterior crown	6° cone secondary*	6° cone secondary*	Telescope secondary*	Telescope secondary*	RS circular*	8-unit bridge
Modellation material:	Wax	Wax	Adapta	Wax	Adapta	Adapta	Pattern	Adapta	Pattern	Pattern	Wax
Alloys: examples in brackets ()											
High-gold content C + B alloys (Degulor M)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
Reduced C + B alloys (Stabilor NF IV)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
High-gold content ceramic bonding alloys (Degudent U)	60–70%	65–75%	65–75%	65–75%	65–75%	50–60%	40–50%	60–70%	55–65%	50–60%	65–75%
Reduced ceramic bonding alloys (Deva 4)	–	65–75%	65–75%	65–75%	65–75%	–	–	–	–	55–65%	70–80%
Palladium-based alloys (Degupal G)	–	70–80%	70–80%	70–80%	70–80%	–	–	–	–	60–70%	75–85%

* We do not recommend rapid heating for preheating of complex restorations.

4 Expansion control

Deguvest Impact is mixed with Deguvest Liquid or with Special Liquid for rapid heating. When these liquids are diluted with distilled water, the expansion can be controlled to compensate for the contraction of the particular casting alloy used.

The higher the concentration of the mixing liquid, the greater the total expansion of the investment.

5 Preparing the casting ring

- Dilute Deguvest Liquid to the desired concentration by adding **distilled water**.
- Line the casting ring with a moistened DeguDent casting ring liner or Deguvest Vlies. The end of the casting ring liner should be flush with the end of the casting ring opposite the sprue base former.
- Position the wax pattern in place with wax on the DeguDent sprue base former.
- Coat the wax model with DeguDent wetting agent Waxit® using a brush or by spraying and remove excess wetting agent (do not use compressed air for this purpose).

- Measure Deguvest powder and liquid in the ratio of 100 g powder : 23 ml mixing liquid resp. according to the following table:

Casting ring	Deguvest Impact powder	Deguvest-mixing liquid
1 x	1 x 60 g = 60 g	14 ml
3 x	1 x 150 g = 150 g	35 ml
6 x	2 x 150 g = 300 g	69 ml
9 x	3 x 150 g = 450 g	103.5 ml

5.1 Mixing the investment material

Premix the powder and liquid for approximately 30 s by hand with a spatula in a clean bowl supplied with the Multivac® mixing unit. Allow the investment to remain under vacuum for about 30 s. Finally, mix under vacuum with the Multivac for about 60 s.

Please use program number 4 for mixing with the Multivac compact (Deguvest Im).

During mixing of the investment, the following points, must be observed:

- High degree of vacuum in the mixing bowl determines the quality. Be sure to wait until the full vacuum is indicated to begin mixing.
- The mixing ratio influences the material characteristics; increasing the amount of liquid has the following effects:
 - thinner consistency
 - extended processing time
 - slightly lower setting expansion
 - slightly lower strength;
 a small amount of liquid added has the following effects:
 - smoother surfaces of the casting
 - somewhat higher expansion;
- Mixing ratio that differs from the recommended range:
 - danger of formation of cracks
 - rougher surface of the casting.
- The mixing time is important for consistent results; strictly adhere to the recommended mixing time.
 - mixing times shorter:
 - uneven casting surface and roughness
 - air bubbles can remain in the investment and result in surface bubbles;
 - longer mixing time:
 - lower setting expansion
 - shorter setting times.

5.2 Investing

Let the mixed investment material flow into the casting ring under light vibration and carefully avoid the formation of bubbles.

Immediately turn off the vibrator after the casting ring has been filled.

Heat treatment	Size of the casting ring			
	1 x	3 x	6 x	9 x
Setting (time)	20 min	30 min	40 min	50 min
Wax elimination (temperature)	300 °C **			
Wax elimination (time)	30 min	40 min	50 min	60 min
Preheating temperature	300 °C **			
Holding time (final temperature) *	20 min	30 min	45 min	60 min

* alloy-specific

** When the casting rings are placed in a preheating furnace preheated to 300 °C, a slow distribution of heat takes place inside the casting ring because of the thermal conductivity of the investment material.

5.3 Setting Time

Let the casting ring set for 20 – 50 min (see table: Heat treatment).

6 Preheating and Casting

6.1 Wax elimination

Place the casting ring in a furnace preheated to about 300 °C and eliminate the wax in about 30–60 min (see table: Heat treatment).

To avoid incomplete cast objects by wax residue within the mould, start the wax burnout immediately after setting of the investment.

6.2 Preheating

Heat the casting ring to the final temperature required for the particular alloy. Maintain the final temperature for a minimum of 20 to 60 min (see table: Heat treatment).

During the production of complex restorations, e.g. large bridges, superstructures or attachments, a heating rate of 7 °C/min is recommended.

The parameters that influence the casting such as wax, resin, adaptic foil etc., shape and size of the restorations, size of the casting ring and type of alloy must be observed individually. All of those factors will influence the final result:

Holding times do not change these factors.

6.3 Rapid heating

Due to its material characteristics, Deguvest Impact is suitable for the rapid heating.

For this purpose, the special liquid must be used. In order to obtain consistent results, the following information is important:

- The use of a **metal casting ring is required**.
- The maximum casting ring size is 6 x and **must not be exceeded**.
- Line the casting ring with a ring liner containing ceramic fibres, e.g. Deguvest Vlies.
The casting ring liner should be flush with the upper rim of the casting ring.
- For high-gold content alloys dilute the special liquid to approx. 50 %.
- For palladium-based alloys dilute the special liquid to approx. 75 %.
- After the investment has set (**after approx. 15 min**), the casting ring should immediately be placed in the furnace which has been preheated to the required temperature (700 – 900 °C). Please ensure that the casting ring is slightly tilted when placed in the furnace.
- In order to achieve homogeneous heating of the casting ring, the holding time at the final temperature is essential (see table: Heat treatment).

Complex restorations should not be heated up rapidly.

6.4 Casting

Cast according to the processing instructions from the manufacturer of the alloy.

6.5 Devesting

After casting, let the casting ring cool down to room temperature and remove the investment.

Soak the investment ring in water for approximately 15 min. Carefully split it open at several points with plaster pliers. When this method is used, the dental casting can be easily removed from the investment, and formation of excess dust is avoided. Final residues of investment can be removed or carefully sandblasted with aluminium oxide, grain size 100 – 150 µm.

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DeguDent GmbH