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Deguvest[®] L

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Deguvest[®] L

- D** Gebrauchsanweisung 3
- GB** Instructions for use 6
- F** Mode d'emploi 8
- I** Istruzioni per l'uso 11
- E** Instrucciones para el uso 13

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Deguvest L is a precision soldering investment developed on a scientific basis. The total expansion of Deguvest L is matched to the dimensional behavior shown by dental precious metal alloys during the soldering process. Deguvest L therefore produces soldered objects with extremely high accuracy of fit.

Deguvest L has outstanding strength. The compound can be shaped well and easily cut, and has high thermal stability. For this reason, Deguvest L can also be used for soldering of metal-ceramic works with high melting alloys.

1 Technical data

Mixing ratio:	depending on the desired consistency, recommended mixing ratio 100 g powder : 20–24 ml water
Total expansion up to 700 °C:	approx. 1.2 % linear
Compressive strength after 2 hours:	3.5 MPa
Processing time span:	3–4 min

2 Instructions for use

- 2.1 Roughen and clean the surfaces of the object to be soldered. The soldering gap should be parallel and should be ≤ 0.2 mm.
- 2.2 Pay attention to correct position of the object on the model. Close the soldering gap with dental wax or modelling wax. Then fix the object with adhesive wax and/or with acrylic which burns without residue. Objects fused with dental porcelain must have no contact with the solder investment compound. Cover porcelain, especially towards the cervical margin and to the basal surface with wax. After the wax has been boiled out, a space is left between the solder investment compound and the porcelain in this way.
- 2.3 Fill Deguvest L and water (distilled water) in the recommended mixing ratio into a clean mixing bowl and thoroughly mix or knead with clean spatula.
- 2.4 Form the solder investment compound into bars and place the object to be soldered with slight hand vibration carefully into the investment compound.
Keep the soldering gap free of investment to avoid sulfur damage with sensitive alloys.

- 2.5** Cut back the soldering block. Remove the adhesive wax or acrylic and boil out the soldering gap. Pay attention when boiling out that no adhesive wax is swept into the soldering gap.
Do not remove the wax with a steam cleaning unit.
- 2.6** Place flux in the soldering gap and on the adjacent metal surfaces.
When objects fused with dental porcelain are soldered, the flux should not touch the porcelain.
- 2.7** Pre-dry the soldering block in the preheating furnace and preheat at approx. 600 °C for 10 min.
- 2.8** Remove the soldering block from the preheating furnace and
- heat up with the soldering torch uniformly to the working temperature of the solder so that the solder flows well into the soldering gap. An excessive overheating of the soldering block must be avoided (danger of sulfur damage).
 - Place the soldering block into the porcelain furnace, apply solder with flux and heat for approx. 1 min without vacuum. Please observe manufacturers' soldering recommendations!
- 2.9** When soldering after the porcelain fusing, and when using alloys with a low coefficient of thermal expansion $CTE (25-600\text{ °C}) \leq 14.2\ \mu\text{m/m}\cdot\text{K}$ (i.e. Degudent® U, Degudent® U94, Deva® 4, Degupal® U and Bond-on 4) the temperature range of 800–900 °C should be avoided. When using all other alloys with a coefficient of thermal expansion $CTE (25-600\text{ °C}) \geq 14.3\ \mu\text{m/m}\cdot\text{K}$ soldering can also be carried out in the temperature range mentioned above.
- 2.10** Let the soldering block cool down to room temperature, remove the object and pickle for approx. 10 min.

3 Safety information

Mixing with other liquids or materials during processing health hazards as well as impairments to the casting quality cannot be excluded.

All dental investment compounds contain parts of crystalline quartz and/or cristobalite.

These substances can adversely affect health if they are inhaled. Can cause silicosis. Always use a sufficiently powerful vacuum air suction or wear a mask.

Carry out the "coarse" removal of the casting after the casting operation under water or soak the investment material mould.

The Safety Data Sheet includes further details.

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