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GA-400-01 AB / GA-400-02 AB / GA-400-02 AC

Operating Instructions

Oil lubricated vacuum pump

GA-400-01 AB / GA-400-02 AB /


GA-400-02 AC

Instructions de Service

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GA-400-01 AB / GA-400-02 AB /

GA-400-02 AC

DENSPLY

Vorwort

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Preface

These operating instructions are valid for oil lubricated rotary vane vacuum pumps type GA-400-01 AB, GA-400-02 AB and GA-400-02 AC.

The operating manual is only valid so long as your equipment complies with the status described within.

These instructions contain all details required for the transport, installation, operation, shutdown and maintenance of these appliances.

Therefore please read the operating instructions carefully prior to the first operation, thus ensuring the safe and economic application of the appliance.

When a fault occurs or maintenance is required which is not dealt with in the operating manual, you should contact our qualified personnel.

All service and maintenance work must be carried out by qualified personnel.

If service and maintenance work is neglected or carried out improperly, our warranty will be invalidated.

Should you have any problems understanding the operating instructions our qualified personnel will be pleased to help you.

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The present operating instructions give you all the necessary advice for the proper and safe handling of the oil lubricated vacuum pumps range GA.

Please read and carefully observe the operating instructions to avoid mistakes and dangerous situations.

Contents

1.	Introduction and general advice	23
1.1	Symbol explanation and definitions	23
1.1.1	Symbol explanation	23
1.1.2	Definitions	23
1.2	General safety advice	24
2.	Technical Data	26
3.	Dimensions	27
4.	Description	28
4.1	Equipment description	28
4.1.1	Agreed application	28
4.2	Function description of the vacuum pump.....	29
5.	Installation of the vacuum pump	30
6.	First operation	30
6.1	Removal of transport safety devices	30
6.2	Fill with oil	30
6.3	Electrical connections	30
6.4	Mount vacuum hose with filter	31
7.	Start-up and shut-down of vacuum pump	31
8.	Motor protection	31
9.	Maintenance	32
9.1	Check oil level.....	32
9.2	Oil lubrication	32
9.2.1	Oil change	32

9.2.2	Execution of oil change	32
9.2.3	Oil changing kit	33
9.3	Care / Cleaning.....	33
9.4	Check the intake filter.....	33
10.	Maintenance.....	33
11.	Trouble shooting and remedy	35
12.	Technical advice and spare parts service	37

List of Illustrations:

2.1	Technical data	26
2.2	Performance diagramme	26
3.1	Dimensions of vacuum pump	27
4.1	Vacuum pump	28
6.1	Connections of vacuum pump	30
6.2	Connection vacuum hose	31
10.2	Disassembl of vacuum pump	34

1. Introduction and general advice

1.1 Symbol explanation and definitions

1.1.1 Symbol explanation

In the manual and on the appliance hand symbols and pictograms are used whose meaning you should memorise. These symbols help you to understand the information in this manual quicker and make you aware of potential danger or important advice.

Attention! Danger sign. It points out danger as well as rules and bans to prevent personal and /or danger to property.

Advice! Attention is brought to advise on the handling and economic use of equipment.

Warning against electric energy! Your life could be in danger. Make sure that all electrical work is carried out by a qualified electrician.

Warning advice. The appliance can start without warning!

Warning of hot surfaces! There is a danger of burns! Observe especially that these surfaces may still be hot after the appliance has been switched off. Work near these surfaces only after they have cooled off.

CE Conformity symbol. It can be recognised by this conformity symbol that the appliance corresponds to the safety guidelines of the European Union.

1.1.2 Definitions

User, Operator: Person who has the authorisation to use the machine. The user must have been given instructions on the safe handling of the appliance from the appropriate authority.

Operator Authority: Responsible for the safe installation, regular servicing and cleaning of the appliance.

Qualified personnel: Trained by the operator authority or by qualified DeguDent personnel who are aware of the dangers of the appliance and familiar with the techniques of the appliance. Qualified personnel is trained to carry out service and repair on the appliance.

1.2 General safety advice

During use, care and maintenance of the appliance, the following fundamental safety measures must be observed for the protection of the operator and qualified personnel as well as the appliance:

 During the development and manufacture of the appliance, the recognised regulations of the technique, as well as the recognised valid standards and guidelines were taken into account and used. In addition the appliance has been designed and constructed in such a way that endangerment through the agreed use are extensively excluded. Nonetheless we feel obliged to describe the following safety measures so the remaining dangers can also be excluded.

Warning! When electrical equipment is being used, the basic safety precautions must be followed, to prevent risk of fire, electric shock and personal injury.

Therefore please read the manual to hand, prior to starting work. Keep this manual within reach for the qualified personnel and the user. The information should be passed on to any successor.

During operation of the appliances the relevant laws and regulations for the place of deployment must be observed! In the interest of safe operation the authorised operator and the supervisory personnel are responsible for keeping within the regulations.

Check during all work on the appliance for possible dangers. All parts must be correctly fitted and all requirements fulfilled to ensure faultless operation. Should the appliance be damaged in any way, the appliance should no longer be used. The appliance should only be repaired by experts. Mark the defect clearly and pull out the mains plug, so that until the repair is effected, no accidents or damage can develop through the defective appliance.

Take account of **environmental influences!** Do not operate the appliance in a wet or damp environment.

Do not use the power connection line for purposes it was not meant for. Do not pull the plug out of the socket by the cable. Only pull the plug from the socket by pulling the plug casing. Protect the cable against heat, oil and sharp edges.

In dangerous situations or during technical problems, separate the appliance at once from the power supply (pull the plug).

Check regularly the power supply line and the casings of electrical components and if they are damaged, have them repaired by a qualified electrician. Check the electrical equipment for external damage before work starts. Check especially whether lines or cable are damaged.

If there is damage, do not operate the appliance. **Pull out mains plug!**

For all maintenance and repair work the mains plug must be pulled out of the socket.

Attention! Work on electrical equipment must only be carried out by a qualified electrician.

Only original spare parts must be used. Otherwise there may be injury to the user.


Warning! The use of other spare parts and other accessories, as stated in this manual, can cause personal injury. Only use spare parts permitted by the manufacture!

2. Technical Data

Type	Suction		Vacuum hpa	Frequency Hz	Fuse (A)	Duration of operation (% ED)	max. ambient temperature (°C)	Sound pressure level db(A)	Weight (kg)	Oil capacity cm ³	Revolution min ⁻¹
	(m/h)	(l/min)									
GA-400-01 AB 230 V	2,4	40	5	50	1,25	100	40	40	10	130	1430
	2,8	48	5	60	1,25	100	40	40	10	130	1700
GA-400-02 AB 100-115 V	2,4	40	5	50	2,15	100	40	43	10	130	1430
	2,8	48	5	60	2,15	100	40	43	10	130	1700
GA-400-02 AC 115-127 V	2,4	40	5	50	2,9	100	40	43	10	130	1430
	2,8	48	5	60	2,9	100	40	43	10	130	1700

fig. 2.1 Technical Data

Vide

fig. 2.2 Performance diagramme

3. Dimensions


fig. 3.1 Dimensions of vacuum pump (mm)

4. Description

4.1 Equipment description

4.1.1 Agreed application

Oil lubricated vacuum pump model GA:

The appliance is intended for aspiration of air and non-aggressive gases. The oil lubricated vacuum pumps are designed for aspirating air from furnaces.


The aspiration of fluids, aggressive or explosive gases is forbidden! It will cause danger to health and the possible danger of explosion or fire!

The appliances are ready-to-connect.

The appliance is designed for operation in a dry, well ventilated room.

The appliance should not be operated in a damp or wet environment. In addition, operating near gases or flammable fluids is forbidden.

- | | |
|---|------------------------------|
| 1 | Silencer
(exhaust filter) |
| 2 | Suction nozzle |
| 3 | electric motor |
| 4 | Pump casing |
| 5 | Pump cover |
| 6 | Oil filler with sealing cap |
| 7 | Oil chamber |
| 8 | Oil sight glass |

fig. 4.1 Vacuum pump

4.2 Function description of the vacuum pump

The oil lubricated vacuum pumps belong to the group of rotary vane vacuum pumps. The pump parts are directly flanged to the electric motor.

The standard built-in winding thermostat prevents overloading of the electric motor. When the winding thermostat is activated, the vacuum pump switches off. The machine will restart automatically, when it has cooled down.

Advice! At unduly high temperatures the appliance will switch itself off.

Warning advice. After cooling off the appliance starts automatically.

The pump part consists of a pump casing (4), a rotor with two vanes, pump cover (5), and an oil chamber (7) with silencer (1).

The brass rotor is positioned off centre in a cylindrical casing. The movable vanes are arranged in lines. During rotation they will be pressed against the wall of the pump casing by centrifugal force. The volume between the vanes alters continuously during rotation. The vacuum is created by this volume alteration.

The well conceived interior construction of the aluminium oil chamber ensures that only a negligibly small amount of oil mist in the exhaust air will reach the surroundings. The oil mist is retained successfully in the big volume oil chamber equipped with cooling ribs and is dripping back into the oil bath.

The pump aggregate is mounted onto a base plate with a carrying handle. In the motor terminal box there is a mains lead with fitted plug. The vacuum pump is switched-in / off by means of the control equipment of the furnace. To reduce vibration the base plate is fitted with four vibration absorbers.

5. Installation of the vacuum pump

Warning! The vacuum pump must be installed lower than the furnace.

6. First operation


fig. 6.1 connections of vacuum pump

6.1 Removal of transport safety devices

- Unscrew socket (2, fig. 6.1) and screw the enclosed sound absorber into thread.
- Remove protective cap (3) from intake socket.

6.2 Fill with oil

- Unscrew cap (1) from oil filling nozzle
- Fill in pump oil up to mark "max" from bottle supplied
- Screw-on cap (1) again.

6.3 Electrical connections

Prior to starting up the vacuum pump check power supply and frequency. Refer to operating instructions for Dentsply ceramic furnaces "First Operation".

6.4 Mount vacuum hose with filter


fig. 6.2 connection vacuum hose

- Push vacuum hose (A, fig. 6.2) onto suction nozzle (3, fig. 6.1) of the vacuum pump.
- Connect hose (E, fig. 6.2) to furnace. Refer to operating instructions for Dentsply ceramic furnaces.

Warning! The non-return valve should be above the vacuum pump.

7. Start-up and shut-down of vacuum pump

The vacuum pump is switched on and off automatically by means of the control system of the furnace.

8. Motor protection

The integrated winding thermostat will switch off automatically, when the maximum permissible ambient temperature is exceeded.

After cooling off, the appliances with temperature switch start again automatically!

If the motor locks for example on starting against pressure or a short circuit in the motor windings, the current consumption can become so high that a fitted temperature switch might get damaged. This can only be controlled in a safe way by the installation of an external overload protector.

When using outside of normal operating conditions the installation of an overcurrent switch is absolutely essential!

9. Maintenance

9.1 Check oil level

The oil level must be kept between marks “min” and “max”.

9.2 Oil lubrication

Warning! Use only recommended pump oil. Do not operate the pump without vacuum for more than 20 minutes (50 hPa).

When using the equipment for short periods the pump will remain cold. This may result in the oil being contaminated with water condensate. In order to ensure reliable operation, comply with the following points:

Continue to run the pump for about an extra 15 minutes, or once a week warm up at minimum vacuum (for at least one hour).

9.2.1 Oil change

The oil change intervals depend on various parameters, such as running time, vacuum level, ambient temperature (min. 12°C, max. 40°C), degree of pollution, water condensate.

Warning! Irrespective of the utilization time, the oil must be changed after every three months at the latest. Otherwise the pump's lifetime will be reduced significantly.

9.2.2 Execution of oil change

- a. Switch-off vacuum pump and remove vacuum hose from nozzle (3, fig. 6.1)
- b. Screw-off sealing cap (1, fig. 6.1) from oil filling nozzle. Remove sealing gasket
- c. Push hose onto oil filling nozzle and hold over oil collecting basin

Warning! Collect the used oil in appropriate container and dispose off under observation of regulations.

- d. Switch-on the vacuum pump; if necessary close air exhaust outlet (4). The used oil is forced out of the oil filling nozzle (1).
- e. Switch-off vacuum pump.

Warning! Refill with oil quantity as stated. Too much lubricating oil will lead to damages. When using not enough oil the lubricating film will be interrupted and this will lead to damages in the pump.

- f. Refill with new lubricating oil up to marking “max”.
- g. Screw on sealing cap.

9.2.3. Oil changing kit

A quick, clean and troublefree oil change is possible with the help of our oil changing kit.

The oil changing kit (article no.: 03107350) comprises:

¼ l special pump oil (for two oil fillings)

¼ l empty bottle to collect the used oil

2 hoses to pump out the used oil (the hose in use is simply put into the empty bottle).

9.3 Care / Cleaning

The vacuum pump should be regularly cleaned on the outside with a cloth soaked in cleaning agent, paying particular attention to the cooling fins, the inspection glass of the oil chamber, and the fan.

If there is hardened lubricating oil or dirt inside the pump, the vacuum pump can be cleaned inside by a cleaning process. Please proceed as follows:

- a. Drain of the oil (see 9.2.2)
- b. Fill in cleaning spirit or petroleum up to the marking “min”.
- c. Run the pump without vacuum, i.e. with the intake open, for about 30 seconds, while closing off the intake opening briefly several times.
- d. Switch off the pump and leave it at a standstill for approx. 30 minutes.
- e. Drain off the cleaning spirit or petroleum and fill in lubricating oil immediately.
- f. Run the pump with full vacuum for approx. 30 minutes.

9.4 Check the intake filter

The transparent intake filter with non-return valve (fig. 6.2). If this filter is soiled, replace the entire filter set not 03260049.

10. Maintenance

Caution! Prior to repair work the vacuum pump must be disconnected from the power supply mains, all lines must be free of pressure.

Repair work must be carried out by qualified specialists.

Warning! The pump casing and the rotor must not be detached from the motor, as they have been aligned in the factory.

Repair work at the parts listed under pos. 1 - 10 may be carried out by service technicians. Other repair works can be carried out only at the factory.


fig. 10.2 Disassembled vacuum pump

Place vacuum pump on working surface and disassemble it following the position numbers (pos. 1 to 10) as shown in fig. 10.2.

Clean all parts using cleaning spirit or petroleum and blow through with compressed air. All parts should have a close inspection. Replace all defective components.

Caution! When working with cleaning agents eye and breathing protectors plus protective clothing must be worn.

Equipment or component parts being sent to the manufacturer for repair, must be carefully wrapped and packed prior to transportation.

Reassembling of the vacuum pump is made in the opposite order. After assembly the exterior surface of the gasket (3, fig. 10.2) has to be sealed with clear lacquer.

11. Trouble shooting and remedy

Warning! Prior to repair work at the pump, plug must be disconnected from power supply!

Failure	Cause	Remedy
<p>Pump does not start</p>	<ul style="list-style-type: none"> • no voltage • pump too cold (<12°C) therefore oil too thick • pump too hot, thermostat switching off • motor fails • pump rotor is blocked, e.g. by dirt in the pump 	<ul style="list-style-type: none"> • Check mains and equipment fuse, if necessary call electrician • warm up pump, use recommended oil • let pump cool down, if necessary allow for better ventilation • replace vacuum pump • have the pump repaired, install suction filter
<p>Pump does not reach performance</p>	<ul style="list-style-type: none"> • Connecting hoses or connections leaking • insufficient oil lubrication and oil seal unsatisfactory • oil contaminated (condensed water, solvent, etc.) • wear in the pump, or shaft seal defective • pump out of action for a long time, oil gums • wrong oil • sinter filter in oil chamber is restricted 	<ul style="list-style-type: none"> • check hoses or connections (seal or replace) • refill oil • change oil, install suction filter • have it repaired • clean pump (see 9.3) • use recommended oil • clean filter (dismantle oil chamber), install new oil chamber seal

Failure	Cause	Remedy
Pump too noisy	<ul style="list-style-type: none"> • not enough oil, lubrication unsatisfactory • pump out of action of a long time, oil gums • bearing damaged • dirt in the pump • ventilator touches ventilator cap • wrong oil 	<ul style="list-style-type: none"> • fill up oil • clean pump (see 9.3) • have motor repaired • clean pump, install suction filter • reposition ventilator cap • use recommended oil
Oil in the suction conduit	<ul style="list-style-type: none"> • operation executed without check valve • check valve leaking 	<ul style="list-style-type: none"> • install check valve in suction conduit • replace it
Oil consumption too high	<ul style="list-style-type: none"> • ambient temperature too high (> 40°C) • wrong oil 	<ul style="list-style-type: none"> • care for better ventilation • use recommended oil
Pump loses oil	<ul style="list-style-type: none"> • oil trickles out of motor (shaft seal defective) • oil chamber surface damaged • oil chamber seal damaged • seal between pump cover and casing, or seal between casing and motor not in order 	<ul style="list-style-type: none"> • have it repaired • replace oil chamber and seal • replace seal • have it repaired

12. Technical advice and spare parts service

Technical advice

DeguDent GmbH

Bürgermeister-Otto-Knapp-Str. 45

49155 Bohmte

Telefon 0 54 71 / 8 07 - 49

Telefax 0 54 71 / 8 07 - 79

Spare parts service

Orders for spare parts should be made quoting spare parts numbers from available parts list and sent to the above address or by telephone:

Telefon 0 54 71 / 8 07 - 49

Telefax 0 54 71 / 8 07 - 79

For ordering spare parts the following details are required:

Model number of vacuum pump and serial number

Position number as in spare parts list

Quantity required

Correct address of the consignee

dispatch details

Terms of delivery

ex works

Repairs

When returning goods, please use the original packing wherever possible.

Always pack the vacuum pump in a plastic bag. Only use recyclable packaging material.

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