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Oxynon®

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Xn

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Oxynton

Oxynton is a special flux for use when soldering non-precious dental alloys.

Non-precious alloys may also be soldered to dental alloys with a high gold content.

Oxynton has very good linking characteristics, ensuring the complete coverage of your working surfaces with solder.

Directions for use:

1. Clean the soldering surface of your object, roughen with a burr or sand blast, whichever is more effective.
2. Fix your objects in the required position using hard wax or modelling resin. Form your soldering block using Deguvest L soldering investment.
3. After melting off the wax with boiling water, the adjoining surfaces should be coated with Oxynton.
4. Heat up the workpiece and apply the solder when the working temperature has been reached. When soldering non-precious and precious alloys together, we recommend using Degulor and Unilot solders. If soldering Biosil to Biosil then use Degudent solder N1W or DWL-Lot 801.
5. Any Oxynton still on the workpiece can be removed by placing the object in a Neacid pickling agent.
6. Always close the lid tightly after use.

Attention!

Harmful if swallowed. Keep out of the reach of children. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable gloves and eye/face protection.

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DeguDent GmbH