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Instructions for use

DENTSPLY

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1 Description

ReproStar is a high-quality addition silicone for the fabrication of duplicate models in the dental laboratory.

This two-component silicone is suitable for all current duplicating techniques and is characterized by:

- Reproduction in great detail
- Stable forms and true reproduction
- High tear resistance
- Thin-flowing consistency
- Stable and smooth cast surfaces
- Simple and safe processing at a 1:1 mixing ratio
- Time saving compared to gel-based duplicating processes

2 Technical specifications

Color component A	White
Color component B	Blue
Final color	Homogenous light blue
Mixing ratio	1:1
Mixing time (manual)	Approx. 45 s
Processing time*	Approx. 5 min
Setting time*	Approx. 30 min
Shore A hardness	14 – 16

* At room temperature (approx. 23 °C). Times may vary at extreme temperatures.

CE 0124

3 Processing notes

3.1 Indications

ReproStar is suitable for the duplicating technique using duplicating flasks.

3.2 Model preparation

The master model and doubling flasks must be **dry** and **clean**. Do not water the model before duplicating.

Note:

Silicone duplicating compounds are susceptible to catalyst poisons such as isolating compounds, cyanoacrylates or resins. These may retard or even prevent setting. Lubricating the duplicating flasks may also result in a fatty layer that does not set.

4 Mixing ReproStar

ReproStar is always mixed at a ratio of 1:1 (component **A**, white : component **B**, blue). If this ratio is not observed, this may adversely affect the setting process and the physical properties of the material.

ReproStar can be mixed using different methods:

- Manual mixing
- Mixing with suitable dosage devices
- Mixing in a vacuum mixing unit

4.1 Manual mixing

Add the ReproStar compounds to a dosage and mixing cup at a ratio of 1:1 and mix well until a homogenous mass is obtained. Pour the compound into the duplicating flask as a thin stream of liquid.

4.2 Mixing with suitable dosage devices

When using 6-kg containers with dosage devices or mixing units, care must be taken to ensure that the minimum filling level is maintained at all times. Otherwise, air will enter the system during processing, which prevents the creation of a homogenous mixture.

4.3 Mixing in a vacuum mixing unit

When mixing the silicone in a vacuum mixing unit, please note the following special instructions:

Use a separate mixing bowl for mixing silicone duplicating compounds. Best results are obtained with an all-new bowl that is free of scratches. Allow excess silicone to harden inside the bowl after mixing. The hardened silicone will be easy to pull away from the edge of the bowl.

Make sure to observe the recommended mixing parameters.

These were defined to protect your unit and to achieve the best possible mixing results. In particular, exceeding maximum recommended **filling levels** may lead to filter clogging and, in extreme cases, to silicone being aspirated into the vacuum pump.

It is recommended not to initiate the mixing process until the requisite vacuum level is reached. Protracted evacuation does not have any effect on the mixing result.

4.4 Setting and deflasking

The silicone will have set after a minimum setting time of 30 minutes (approx. 23 °C). At that time, you can carefully remove (deflask) the master cast from the silicone mould.

CAUTION:

Removing the master model will distort the silicone mould. However, the silicone is capable of elastic recovery and will return to its proper shape within a few minutes. It is recommended not to initiate fabrication of the refractory model until the silicone has rested for 15 minutes.

4.5 Pouring the silicone mould

The silicone mould is highly elastic and permanently retains its exact shape. Before pouring with investment compound, the silicone mould can be wetted with a surface-tension relieving agent to improve the flow of the investment compound on the silicone surface. Carefully blow-dry the mould. Check that the silicone surface is dry and completely set.

Beside of investments for pouring out all dental plasters and model-acrylics can be used.

4.6 Disinfection

ReproStar should be stored at room temperature; only then can the minimum shelf life of 2 years from date of production be ensured.

5 Storage conditions

Store ReproStar at a temperature between 15 °C (59 °F) and 28 °C (82 °F).

6 Presentation

ReproStar component **A**, white
and component **B**, blue

in two 860 ml (1,000 g) bottles 2511 1502

in two 5,170 ml canisters (6,000 g) 2511 1512

www.dentsply-degudent.com

DeguDent GmbH

A Dentsply International Company
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
www.dentsply-degudent.com

DENTSPLY