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F3 | ERGO

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Gebruiksaanwijzing

F3 | ERGO

D	Gebrauchsanweisung	3
GB	Instructions for use	11
F	Mode d'emploi	19
I	Istruzioni per l'uso	27
E	Instrucciones para el uso	35
NL	Gebruiksaanwijzing	43

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Seite 10

Vakat !

Milling machine for producing high precision dental restorations

1 Packing list

1 Milling machine, complete	1 Dust cover
1 Milling spindle	1 Velcro for securing the (milling) arm
1 Power cable	1 Safety sticker
1 Compressed air hose	1 Instructions for use / Guarantee card

2 Technical data

Power supply voltage:	Refer to rating plate
Mains frequency:	50 Hz / 60 Hz
Power input:	150 W
Safety fuse	
for mains voltages 100 – 115 V:	4 A
for mains voltages 230 – 240 V:	2 A
Compressed air supply:	min. 4 bars, max. 8 bars
Speed of milling spindle:	
Clockwise:	1,000 – 25,000 r.p.m.
Anti-clockwise:	1,000 – 7,500 r.p.m.
Magnetic model platform:	Ø 150 mm
Vertical travel:	225 mm
Dimensions (H/W/D):	600 / 270 / 490 mm
Weight:	approx. 25 kg
Operating conditions:	+10 °C – +35 °C
	Must only be used in dry areas
Storage conditions	-15 °C – +35 °C
	5 % – 80 % relative humidity, non-condensing

3 Accessories

Item	REF
Illumination device	5153 0100
Model platform	2351 1033
Multi-purpose spindle	5153 0101
Spacer plate for adjusting the height	5153 0103
Plastering platform with 44 mm Ø rubber sleeve	2351 1035
Plastering platform with 65 mm Ø rubber sleeve	2351 1036
Plastering platform with 80 mm Ø rubber sleeve	2351 1037
Transfer platform for restorations with conical crowns	5153 0102

4 For your safety

Read these operating instructions carefully before installing and operating the milling machine.

Always heed these instructions.

The rotary instrument flings off shavings during milling and drilling. Always put on safety glasses before operating this milling machine.

Always comply with the safety precautions for rotary instruments. Risk of injury.

Only use original spare parts and accessories for the **F3** | ERGO milling machine. This is the only way of ensuring that the performance figures and safety standards are attained.

Do not open the milling machine. There are dangerously high voltages inside this milling machine.

Maintenance work and repairs must only be carried out by authorized service technicians.

To prevent short-circuiting and damage to the milling machine, protect it against moisture.

If the milling machine is not used for a longer period of time, it should be unplugged from the mains supply.

5 Description of the product

The **F3** | ERGO is an electric milling machine for fabricating high precision dental restorations.

6 Installation

The **F3** | ERGO milling machine is a bench unit.

- **Place the milling machine on a sturdy bench to ensure that it can be used safely and precisely.**
- Push the drawer into the base of the milling machine.
- Insert the milling spindle into the holder on the milling arm. **Ensure that the positioning pin, which is visible in the spindle holder, locks into the groove in the milling spindle.**
- Turn the locking nut to lock the milling spindle (A) (Figure 1).
- Plug the milling spindle cable into the socket in the column and tighten the ring on the plug (B) (Figure 1).
- Use the quick-link coupling (C) to connect the compressed air hose to the right side of the milling machine. **Insert the hose to the stop (approx. 10 mm)** (Figure 2).
- Connect the milling machine to your compressed air supply. You will need a compressed air coupling with a nominal diameter of 7.

Caution! The hose must always be connected to the milling machine first as, otherwise, the compressed air would escape from the other end.

- **Before connecting the **F3** | ERGO milling machine**, ensure that the voltage quoted on its rating plate coincides with your mains voltage.

- Plug the power cable into the socket (D) on the right side of the milling machine and connect it to the mains supply (Figure 2).
- Stick the “Wear safety glasses” sticker supplied with the milling machine to an appropriate place on the bench.

The **F3** I ERGO milling machine is now ready for use.

7 Operating instructions

7.1 Electric controls

The **F3** I ERGO is switched on with the power switch (E) on its right hand side (Figure 2). The LED (F) lights up to show the direction of rotation (Figure 3).

7.2 Milling spindle – Chuck

Do not change the rotary instrument until the milling spindle has been switched off and come to a standstill. Opening the chuck clamps the milling spindle and may damage it.

To open the chuck, pull the clamping lever forwards, over the resistance, to the stop. The chuck is closed in the reverse order.

The chuck has not been closed properly until the two spots on the milling spindle housing are aligned.

Fig. 1


Fig. 2



Fig. 3

To replace the chuck, remove the milling spindle from its holder and withdraw the rotary instrument from the chuck. Insert the shank of the instrument into the chuck to prevent it being pressed together. With the chuck **open**, use the chuck key, supplied with the milling machine, to unscrew it from the spindle. Use a second chuck key to prevent the shaft rotating. Fit the new chuck in the reverse order, ensuring that it is tightened properly.

With regard to the service life of the chuck, ensure that an instrument is always fitted to it and tightened – **even when it is not in use.**

Press the  touch controller to switch the milling spindle on. The instrument fitted to the spindle will then rotate at the speed selected on the speed controller. The speed can be adjusted progressively while the instrument is running.

The direction of rotation can also be changed with the  touch controller while the instrument is running. When the instrument is running anti-clockwise, the speed is automatically limited to 7,500 r.p.m.

7.3 Magnetic model platform

The magnetic model platform is switched on and off with the  touch controller.

7.4 Hinged arm

The hinged arm can be adjusted vertically with the  and  touch controllers.

The feed motor switches off automatically once the highest or lowest positions are reached.

The hinged arm is fixed in the horizontal position by pressing the  touch controller to activate the pneumatic brake.

The **F3** | ERGO includes an instrument change mode which is activated with the  touch controller. The instrument change mode raises the hinged arm 50 mm to allow adequate space for replacing the instrument in the milling spindle. The same touch controller is then pressed again to return the hinged arm to its working position. **For safety reasons, the instrument change touch controller must be pressed continually while the arm is moving downwards.**

While the hinged arm is in the instrument change mode, LED (G) flashes and the vertical feed touch controllers are deactivated (Figure 3).

7.5 Drilling unit

Once the locking screw (H) on the left side of the forward feed jig has been loosened, the jig can be lowered with the lever on the right side. The drilling depth can be read off the scale (I) on the hinged arm. The depth stop of the jig is adjusted with the micrometer screw (J) (Figure 4).

The sleeve (L) is pulled toward the knob to unlock the forward feed lever (K) and fix it in its vertical position.

The forward feed jig is sprung vertically. The force of the relieving spring is increased or decreased by turning the threaded rod with the knob. It can be fixed at various levels with the knurled nut (M) (Figure 5).

Fig. 4



Fig. 5

8 Care and maintenance

8.1 Maintenance

The **F3 I** ERGO requires **virtually no maintenance**. The **F3 I** ERGO must only be repaired by **Degussa Dental** employees or a workshop authorized by **Degussa Dental**.

8.2 Care


The **F3 I** ERGO can be cleaned with a soft cloth but without using solvents. Avoid blasting the milling machine with compressed air as this may damage it.

8.3 Replacing the safety fuse

Unplug the milling machine from the mains !

The safety fuse (N) is on the right side of the milling machine beside the power switch (E).

Should the fuse be triggered off several times, the milling machine is probably defective and must be repaired by an authorized service technician.

9 Guarantee terms

Every milling machine is tested for perfect functioning prior to delivery.

We guarantee this milling machine under the following terms:

The guarantee period begins on the date of delivery and expires 12 months later. Any defects due to faulty materials or manufacturing, which occur during the guarantee period, will be repaired free-of-charge. No further claims can be taken into account.

Carriage and packaging charges as well as customs and import duties in the purchaser's own country will be charged to him unless the law demands otherwise.

The guarantee is null and void if repairs or other interventions are carried out by persons not authorized to do so by us. Claiming against the guarantee does not prolong the guarantee period.

Guarantee claims can only be accepted if the guarantee certificate is presented **together** with the invoice or delivery note.

For use in dentistry only.

10 Declaration of conformity

EU Declaration of Conformity, EU Machines Directive (98/37/EG)

Degussa Dental

Degussa Dental GmbH

Postfach 1364 · D-63403 Hanau

We hereby declare that the:

F3 | ERGO milling machine

conforms with the basic requirements of the following EU directives:


1. **EU Machinery Directive 98/37/EC**
2. **EEC Directive on electrical equipment designed for use within certain voltage limits 73/23/EEC**
3. **EEC Directive relating to electromagnetic compatibility 89/336**

The following harmonized standards were employed:

EN 292-1:1991, EN 292-2:1995, EN 60204-1:1998
EN 61010-1:1994, EN 50081-1:1993, EN 50082-1:1997

Torsten Schwafert
Director business area
dental laboratory products

Dr. Alexander Völcker
Research & Development

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