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Spotlight LED

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones de uso

1. Auflage Spotlight LED / REF 5153 0110

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1st Edition Spotlight LED / REF 5153 0110

This documentation is valid for Spotlight LED models with serial numbers from xx00101xx onwards. You can work out the serial number of your device on the nameplate by crossing out the first two and last two digits on the attached numerical code.

Before connecting and using the Spotlight LED, please read these Instructions for Use carefully!

Target group:

This document is intended for users of the milling machines F3 ergo and F4 basic.

Documentation contents:

In this manual you will find information you need for:

- installing and starting up
- operation and handling

Manufacturer:

DeguDent GmbH

Postfach 13 64

D-63403 Hanau

Telephone +49 (0)61 81 59 50

Fax +49 (0)61 81 59 58 58

www.degudent.com

© May 2004 by DeguDent GmbH

We reserve the right to make changes and correct errors. This document may not be duplicated without our permission.

Table of Contents

1	Description of Product	10
	1.1 Scope of Delivery	10
	1.2 Description of Device / Technical Features	10
2	Safety	10
	2.1 Intended Use	10
	2.2 Design in Accordance with EC Directives	11
	2.3 Requirements for the Place of Installation	11
	2.4 Safety Tips	11
3	Getting started	12
	3.1 Connecting the Device	12
4	Care and Maintenance	12
	4.1 Positions for Maintenance	12
	4.2 Care	12
	Declaration of Conformity	13

1 Description of Product

1.1 Scope of Delivery

The Spotlight LED (REF 5153 0110) is delivered with the following:

- 1 illumination unit Spotlight LED
- 2 cross-head screws M3 x 35
- 1 instructions manual / guarantee card

Please keep the transport packaging

The supplied packaging protects the device from transport damage.

If possible, please keep the packaging materials in case you need to transport the device at a later date.

1.2 Description of Device / Technical Features

The Spotlight LED is an illumination unit to be operated exclusively with the F3 ergo and F4 basic milling machines. The Spotlight LED serves to illuminate the working field for the F3 ergo und F4 basic milling machines.

Technical data

Mains voltage	100–240 V
Mains frequency	50/60 Hz
Wattage	3 W

2 Safety

Please read this chapter before starting your device for the first time!

Observe the information contained in these Instructions for Use under all circumstances!

2.1 Intended Use

Only operate the illumination device in its assembled state in conjunction with the F3 ergo und F4 basic milling machines.

Never open the device. Hazardous voltages exist inside the device. Maintenance work and repairs may therefore only be undertaken by authorized service technicians.

Protect the device from penetration of moisture in order to avoid short circuits and damage to the device.

Please disconnect the mains plug if you do not intend to use the device for a long period.

2.2 Design in Accordance with EC Directives

The Spotlight LED complies with the EC Directive 89/336/EEG (Electromagnetic Compatibility).
The following standards were used to comply with this directive:

- EN 55015: 2000 + A1: 2001 + A2: 2002
- DIN EN 61000-3-2: 2000

The following standards of EC Directive 73/23/EEC (Low Voltage Directive) were used:

- EN 61347-1: 2001
- EN 61347-2-2: 2002
- EN 61547: 1995+A1: 2000
- EN 60598-1+ A12: Nov. 2002

2.3 Requirements for the Place of Installation

The device may only be operated continuously in dry rooms.

It may not be operated continuously outdoors.

This device may not be operated in an atmosphere where there is a risk of explosion.

It may not be operated when there is a lot of dust (or even abrasive dust).

Ambient temperature permitted for operation:

Allowable storage temperature range: -15 °C – +40 °C

Allowable relative humidity (without condensation):

max. 90 % at +40° C

Connecting to the mains

The Spotlight LED can only be connected to the milling machines F3 ergo or F4 basic in its assembled state.

Under no circumstances should you attempt to make modifications to the machine or the power supply cable supplied with it.

Your power supply connection must have a functional protective earthing.

2.4 Safety Tips

Only an authorized service technician may perform repairs on the device

Repair work not carried out properly can impair the safety and protective functions (including electro-magnetic behaviour) of the Spotlight LED, which may be hazardous to health or damage the device.


3 Getting started

3.1 Connecting the Device

Please follow the procedure below:

Installation

- Remove the mains plug on the milling machine.
- Fix the casing (1) to the rear of the milling machine with the two screws (2) supplied as shown in the illustration.
- Join the milling machine and illumination unit with the bridging cable (3).
- Connect the illumination unit with the mains cable (4).

3.2 Operating

- Switch on the lamp by activating the toggle switch.
- Adjust the flexible light holder as required.
Do not go below the minimum bending radius of 50 mm.

4 Care and Maintenance

4.1 Positions for Maintenance

The Spotlight LED requires no maintenance. Only DeguDent GmbH staff or a workshop authorized by DeguDent GmbH may carry out repairs to the Spotlight LED.

4.2 Care

The Spotlight LED can be cleaned with a soft cloth or by the use of solvents. Please do not air-abrade the device, as this may damage it.

Declaration of Conformity

DeguDent GmbH

Rodenbacher Chaussee 4
D-63457 Hanau-Wolfgang

We hereby declare that the:

Spotlight LED

conforms with the basic requirements of the following EU directives:

- 1. EEC Directive on electrical equipment designed for use within certain voltage limits (Low Voltage Directive) 73/23/EEC**
- 2. EEC Directive relating to electromagnetic compatibility 89/336/EEC**

The following harmonized standards were employed:

EN 61347-1: 2001

EN 61000-3-2: 2000


EN 61347-2-2: 2001

EN 61547: 1995 + A1: 2000


EN 55015: 2000 + A1: 2001 + A2: 2002

EN 60598-1 + A12: Nov. 2002

Hanau, 10. 2. 2004

Dr. Udo Schusser
Director of Research & Development

Torsten Schwafert
Director of Operative Area

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