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# Neacid<sup>®</sup>-Beizgerät

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso

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Vakat !

# Neacid<sup>®</sup>-Beizgerät

<b>D</b>	Gebrauchsanweisung . . . . .	4
<b>GB</b>	Instructions for use . . . . .	7
<b>F</b>	Mode d'emploi . . . . .	10
<b>I</b>	Istruzioni per l'uso . . . . .	13
<b>E</b>	Instrucciones para el uso . . . . .	16

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## 1 Intended use

The Neacid pickling unit was developed for material-specific pickling (deoxidizing) of precious metal work pieces together with the pickling agent Neacid<sup>®</sup> (information on Neacid can be found under section 3 of the operating instructions). The thermostatic control of the electrically heated unit and the resulting low power consumption allow permanent heating during the entire working day. Therefore pickling processes can be carried out with the unit at any time and in rapid succession.

The swivelling sieve that is connected to the unit housing allows to avoid the use of steel tweezers or similar tools for the transport of the work pieces.

The unit is designed for indoor use only. The room temperature should be between 5 °C and 40 °C. Maximum relative air humidity of 80 % for temperatures up to 31 °C, linearly decreasing to 50 % of relative air humidity at 40 °C. External power supply. Voltage variations not larger than  $\pm 10$  % of the rated value. Contamination degree 1 or 2 according to IEC 664 (section 3.7.3).

### **Please note!**

The Neacid Pickling Unit was developed for use with Neacid pickling agent. Other pickling agents may decompose the plastic components, leak into the unit or power switch and cause a short-circuit!

Neacid (Pack of 2 sachets 2 x 65 g)

REF 53 2522 0001

**Please adhere to the operating instructions since operation contrary to the operating instructions may result in hazards for the operator.**

## 2 Technical data

Voltage:	200/240 V~, 50/60 Hz
(special version 100/120 V~, 50/60 Hz)	
Power consumption:	40/60 W at 100/120 V~ 40/60 W at 200/240 V~
double insulation	
Housing, bath container, cover and sieve are made of acid-resistant and shock-proof polypropylene.	
Dimensions:	L 180 x W 170 x H 190 mm
Weight:	650 g
REF	53 5431 0001 (200/240 V) 53 5431 0017 (100/120 V)

## 3 Information on Neacid

Neacid is a pickling agent in powder form. When it is dissolved in water a highly efficient acid bath is formed to remove oxides and fluxes from cast and soldered objects produced from precious metal alloys.

It is also suitable for 14 and 18 carat colored and white gold as well as for silver alloys. **For Palliag® and silver alloys a pickling time of 5min must not be exceeded.** These alloys will not become entirely bright; copper oxides, however, will be dissolved without any residues. During extended pickling these alloys will be attacked and obtain a darker color. For 8 carat decorative gold the risk of stress corrosion cracking exists so that Neacid can only be recommended for stress-free casting.

### Attention during the handling of Neacid!

Eyes, mucous membranes and open wounds must neither get in contact with Neacid solution nor with Neacid powder.

### Irritant

Irritating to eyes and skin.

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Keep out of reach of children.

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

After contact with skin, wash immediately with plenty of water.

Avoid release to the environment. Refer to special instructions / Safety data sheet.

The safety data sheet includes further details.

## 4 Putting into operation

In order to avoid tilting and thus spilling of the pickling liquid, the unit must be placed on a solid and flat surface.

Prior to initial operation fill the bath container with 250 cm<sup>3</sup> of water, connect the unit to the mains supply and switch it on. The control light in the push-button switch lights up.

After approx. 20 min pour a bag with 65 g of Neacid into the heated water and stir with a glass or wooden stick until the powder is completely dissolved. After approx. 30 min after switching on the unit the pickling liquid will have the optimum working temperature of approx. 70 °C.

## 5 Working instructions

**5.1** For pickling lift the cover vertically, place the work piece into the swivelling sieve and close the bath container again.

In order to avoid damage to the swivelling sieve do not insert any hot objects.

Additionally, hot objects may cause splashing of the acid.

**5.2** The pickled work pieces can be rinsed directly at the working place in a container that has been filled with water. During this process, however, the rinsing water must be exchanged repeatedly. Therefore we recommend to rinse the pickled work pieces under running water. This process can be carried out quickly with the Neacid pickling unit. The bath container with the acid is left at the working place and the lower side of the cover is used as a dropping tray during the transport of the work pieces.

**Metal tweezers must not be immersed into the acid** since this may result in discoloration of the casting objects.

After rinsing, take the work piece out of the swivelling sieve and immediately place the cover on the unit again in order to avoid unnecessary evaporation of water from the pickling liquid. Damage to the unit can not be excluded if the unit is continuously operated with a filling level that is too low.

**5.3** If the liquid level has dropped below the marking due to the evaporation of water, the bath must be filled up with water.

The pickling liquid should be exchanged when it displays intensive blue-green discoloration.

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**DeguDent GmbH**