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# Waxtip

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso

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## Brief Description

The Waxtip is a unit designed to heat the modelling instruments known from the Sensorwaxer; it is connected directly to an electric power socket, and is very easy to operate. At the Waxtip – depending on the phase of work – a temperature of 60 °C – 230 °C can be set.

By means of a thin cable, the Waxtip is connected with the particular working instrument attached to it. Both the sensor instruments AS 2 and BS 2, as well as the tools C 2, D 2, and E 2 which operate without a sensor can be used.

## A. Equipment Supplied

1 Waxtip (control unit)

REF 2336 0001

### Accessories

Designation	Design	REF
Instrument AS 2	Waxing-up probe with sensor	2351 1160
Instrument BS 2	Small modelling knife with sensor	2351 1161
Instrument C 2	Large modelling knife without sensor	2351 1162
Instrument D 2	Large wax spoon without sensor	2351 1163
Instrument E 2	Wax preheating ring without sensor	2351 1164

## B. Technical Data

Power voltage/Power frequency:	100–120 VAC 50/60 Hz 200–240 VAC 50/60 Hz 127 VAC 60 Hz
Power consumption:	12 W
Unit fusing:	T 0.16 A
Height/width/depth:	120/70/60 mm
Weight:	0.5 kg

## C. Instructions for Operation

1. Plug the Waxtip into a power socket near the sorking area. The Waxtip is now ready for operation.

### 2. Connection of a sensor instrument (AS 2 or BS 2)

2.1 Insert the diode plug attached to the sensor instrument into the socket on the Waxtip.

2.2 Operate the on-off ("Netz") switch; the green indicator lamp in the switch lights up.

2.3 At the same moment the sensor instrument is heated up to 60 °C.

2.4 Set the working temperature at the temperature controller (from one mark to the next ~ 17 °C).

2.5 When the two sensor rings at the front end of the black handle are now touched **simultaneously**, the working tip is heated up in just a few seconds to the temperature set at the controller; this process is indicated by the red diode lighting up.

2.6 When the two sensor rings are released (red diode switches off), the working tip cools down to the so-called stand-by temperature of 60 °C. This process can be repeated as often as desired.

### 3. Connection of an instrument without (C 2, D 2 or E 2)

3.1 Insert the diode plug attached to the instrument into the socket on the Waxtip.

3.2 Set the desired working temperature between 60 °C and 230 °C (from one mark to the next ~17 °C).

3.3 Operate the on-off ("Netz") switch; the green indicator lamp in the switch and the red diode light up.

3.4 The working tip of the instrument is automatically heated to the temperature that has been set, and it is then held constant at this temperature. The automatic control system that otherwise is used with the sensor instruments remains inactive in this mode of operation.

## D. Special Information

1. The Waxtip is designed for continuous operation; i.e., the unit can remain switched on all day. At the end of a working day, however, the Waxtip should be switched off for safety reasons.
2. The working tips on the individual instruments are produced from a high-grade silver alloy because of the good thermal conductivity. In addition, the ductility of the alloy permits certain corrections in form at the working tips through bending, grinding, or similar processes if this should become necessary for a special intended application.
3. Avoid damage to the cables caused, for ex., by hot tips on the working instruments!

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