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Instructions for Use Solaris Pre-Liquid

(GB)

Solaris Pre-Liquid is a degreasing and activation liquid used for the pre-treatment of prosthetic components to be gold-plated. Solaris Pre-Liquid ensures optimal degreasing and activation of the metallic surfaces.

Amount of liquid 1000 ml Working temperature .. approx. 50-60°C (122-140 °F) Liquid movement Ultrasound
 Duration of the treatment approx. 30 min Colour of the liquid milky Stress relief test water is not repelled

First-aid measures: See EU Material Safety Data Sheet (MSDS)

Protective measures: See EU Material Safety Data Sheet (MSDS)

Dosage: Use 500 ml of liquid in a large beaker. If required, add more Solaris Pre-Liquid until all objects to be treated are covered.

Processing:

Caution: Solaris Pre-Liquid must never be mixed with Solaris Flash-Liquid or Solaris gold baths. Avoid contamination, and thoroughly rinse the beakers under running water after use. For additional processing instructions, consult the "Solaris - Instructions for Use for Friction Renewal and Gold-Plating". Solaris Pre-Liquid is suitable for multiple use. It can be used between 5 and 10 times or until it is visibly contaminated or the degreasing effect is diminished.
 To degreased and activated the prosthetic components, treated them with Solaris Pre-Liquid in an ultrasonic bath at temperatures between 50 and 60 °C (122-140 °F) for approximately 30 minutes. Rinse under running water. The desired degreasing effect has been achieved if the water is not repelled and the water film does not separate. If that is not the case, the degreasing process with Solaris Pre-Liquid must be repeated.

Storage: Return Solaris to the original container after use and close the container tightly.

Use before date: See label

Disposal Any product residue must be disposed of in accordance with applicable laws and regulations.

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DeguDent GmbH