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# CoCr- Modellguss- legierung

REF 3502 0001

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso  
Instruções de utilização  
Instrukcja stosowania  
Kullanmakilavuzu  
Инструкция по  
применению  
تعليمات التشغيل

CE 0124

**DeguDent**

A Dentsply International Company

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D

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# CoCr Model Casting Alloy

This high-quality casting alloy from DeguDent for model casting is characterized by harmonized technical qualities. It facilitates casting of fine structures.

It can be melted down using all high-performance dental material melting systems and processed using tested phosphate-bound investments.

## 1 Contraindications

Do not use in case of hypersensitivity to one or more of the metal constituents of the alloy.

## 2 Side effects

Potential side effects include allergies against metals contained in the alloy and discomforts of an electro-chemical genesis. Systemic side effects due to alloy metals are claimed in isolated cases.

## 3 Interactions

Avoid occlusal and approximal contact between different alloy types.

## 4 Technical data (mean values)

Melting interval	1320 – 1380 °C
Vicker's hardness	400 HV10
Modulus of elasticity	220 GPa
0,2%-Expan. limit	650 MPa
Tensile strength	850 MPa
Rupt. limit	5 %
Density	8.4 g/cm <sup>3</sup>

## 5 Chemical composition (in % by mass):

64.8 % cobalt, 28 % chromium, 5 % molybdenum, silicon, manganese, carbon, nitrogen, nickel, < 0,1 % each.

## 6 Processing information

### 6.1 Casting sprues

Each cast object should, depending on the frame dimensions, have 2 – 4 casting channels: use 3 – 4 mm dia. wax wires at the thickest point, not tapered, to facilitate inflow and controlled solidification.

### 6.2 Investing

Apply thin coat of wetting agent (e.g. Waxit®) to wax object with a brush or spray, then blow dry carefully (do not use compressed air). Invest model with model cast investment compound (e.g. Optivest® or Biosint®-Supra). See information on use of the investment compound. Use of a fine investment compound is possible, but it is not necessary if a fine-grained investment compound such as Optivest® is used.

### 6.3 Pre-heating

Use manufacturer's recommendations for preheating times and temperatures. After the final temperature of 1,000 °C is reached pre-heat the casting cuvette for 60 minutes. For particularly thin constructions, raise pre-heating temperature to 1,050 °C. If larger numbers of casting rings are in the oven, the temperature holding time must be increased accordingly.

#### 6.4 Casting

The CoCr alloy can be cast using all sufficiently high-performance melting and casting equipment.

- a) Melt with the reducing flame of a propane-oxygen or natural gas-oxygen or acetylene-oxygen mixture. Use ceramic melting boxes. Cast after liquidus temperature is exceeded, before the oxide skin ruptures.
- b) When using the HF casting equipment (e.g. Degutron®) do not use a graphite or carbon crucible. Cast before the oxide skin ruptures.

#### 6.5 Disembedding and cleaning of the casts

Following casting, allow to cool for at least 20 minutes. Do not quench in water. Blast with corundum or similar blasting material. To preserve the passivation layer on the surface, do not pickle CoCr alloys.

#### 6.6 Finishing

Finish processing of CoCr model casting alloy is easily done using standard finishing instruments for model cast prostheses. For electrolytic polish, protect braces with cover lacquer from wear.

#### 6.7 Soldering / brazing

Solder connections of CoCr model casting alloy to noble metal alloys using flux (e.g. Oxynon®). The solder to be used depends on the noble metal alloy to be soldered. For repair and connection solderings of the CoCr model cast alloy, we recommend Degudent Solder N1W.

#### 6.8 Re-use

Highly pure materials are required for production of a high quality dental prosthesis. For this reason, only original casting cubes should be used to avoid risks.

## 7 Supply forms

CoCr Model Casting Alloy (1.000 g)  
REF 3502 0001

For dental use only.

#### **Warning!**

CoCr Model Casting Alloy contains cobalt and chromium. There are maximum allowable concentrations (MAC values) for certain compounds containing these substances. Elaboration with rotating instruments should fundamentally be undertaken with object exhauster. Furthermore, we recommend wearing mouth guard or face shield.

In cases of a known hypersensitivity to one of the components, this alloy should not be integrated or machined.

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