# BiOcclus<sup>®</sup> inlay

REF  
1206 0002

50534402/b  
Last revision: Aug. 06

Au and Pt-metals	Au	Pt	Ag	Zn	Ir	
	83.3	4.8	10.0	1.8	0.1	mass-%

DeguDent

GB

High-gold alloy for the fabrication of dental restorations (crown-and-bridge)

Type hard Color yellow

Contraindications: Do not use in known hypersensitivity to one or several metals contained in the alloy.

Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

Interactions: Avoid occlusal and approximal contacts of different alloy types.

Avoid inhaling dust and vapours while in contact with dental alloys. Use suitable vacuum appliances for protection against dust and vapours. Use a facemask or respiratory protection.

For additional information concerning handling, dust protection and application notes, please consult our brochure "Precious-metal dental alloys – Processing instructions".

For dental use only.

Melting range	Preheating temperature	Casting temperature	Annealing	Vickers hardness	0,2% Yield strength*	Tensile strength*	Elongation*	Density
°C	°C	°C	°C min not possible	HV 5 w) 100 g) 145	N/mm <sup>2</sup> w) 230 g) 340	N/mm <sup>2</sup> w) 370 g) 440	% w) 13 g) 4	g/cm <sup>3</sup> 17.4
980-1060	700	1200						

w: soft (15 min/750 °C)

g: Achievable values after slow cooling inside the mould

\* Measurements performed on samples pursuant to DIN EN ISO 1562

Instructions for Use

Modelling: When modelling and finishing the frameworks, make sure to maintain a minimum wall thickness of 0.3 mm for single crowns and 0.5 mm for abutment crowns.

Spruing: Fit the waxed-up bridge framework with adequately dimensioned sprues.

Bar casting: The diameter of the head canal is 3-4 mm. The manifold canal has a diameter of 4-5 mm depending on the volume of the casting. The volume of the bar should correspond to that of the casting. The distance between the manifold canal and the waxed-up framework should be 5-8 mm, diameter 3 mm.

Investing: In phosphate-bonded / gypsum-bonded investment compounds (such as Deguvest<sup>®</sup> SR, F, CF and California).

For rapid heating we recommend Deguvest<sup>®</sup> Impuls, SR and California.

Wax elimination/preheating: At 300 °C and 700 °C depending on flask size. Follow the manufacturer's recommendations.

Casting: The alloy can be safely processed in resistance-heated furnaces in graphite crucibles, or in open-flame casting units in the ceramic-melting trough.

Casting temperature 1200 °C

At least one-third of new material must be used. Sandblast sprues/cones before reuse to thoroughly

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Temperatura di fusione 1200 °C.

Casting: The alloy can be safely processed in graphite-heated furnaces in graphite crucibles, or in open-flame casting units in the ceramic-melting trough.  
 Casting temperature 1200 °C  
 At least one-third of new material must be used. Sandblast sprues/cones before reuse to thoroughly remove all oxides or residual investment.  
 Post-heating: Resistance-heated: 30 s, Prestomat: 15 s, oxygen/propane; HF and arc-casting machine: 5-10 s.  
 Divesting: Allow the flask to cool to room temperature after casting and divest the castings using 50-µm glass beads at no more than 2 bars of pressure away from the investment.  
 Do not use a hammer for divesting!  
 Finishing: Rotating instruments made of hard alloys, ceramic or cemented grinding heads, polishing cloth and rubber polishing tools can be used to process alloys that cannot be fired on.  
 Pickling: After casting or soldering, pickle the objects for a maximum of 15 min in a Neacid® bath.  
 When selecting the solder, please note the temperature interval between the soldering temperature and the alloy's solidus temperature.  
 We recommend Context as an antflux.

Soldering: Lot DG 750 (750 °C)  
 Flux: Anoxan®  
 Laser welding: No laser welding wire is available.  
 Annealing: not possible

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