# BiOcclus® N

| Au and Pt-metals | Au   | Pt   | Zn  | Nb  |        |
|------------------|------|------|-----|-----|--------|
| 97.8             | 85.6 | 11.9 | 1.6 | 0.4 | mass-% |

| Au and Pt-metals | Fe  | Ir  | Rh  |        |
|------------------|-----|-----|-----|--------|
| 97.8             | 0.2 | 0.1 | 0.2 | mass-% |

REF  
1316 0002

50534404/a  
Last revision: Nov. 05

DeguDent

GB

High-gold alloy for the fabrication of dental restorations (for ceramic bonding)  
 Type extra hard Color yellow  
 Contraindications: Do not use in known hypersensitivity to one or several metals contained in the alloy.  
 Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.  
 Interactions: Avoid occlusal and approximal contacts of different alloy types.  
 Avoid inhaling dust and vapours while in contact with dental alloys. Use suitable vacuum appliances for protection against dust and vapours. Use a facemask or respiratory protection.  
 For additional information concerning handling, dust protection and application notes, please consult our brochure "Precious-metal dental alloys – Processing instructions".

For dental use only.

| Melting range | Preheating temperature | Casting temperature | Mean linear CTE |             | Vickers hardness         | 0,2% Yield strength*                  | Tensile strength*                     | Elongation*       | Density                   |
|---------------|------------------------|---------------------|-----------------|-------------|--------------------------|---------------------------------------|---------------------------------------|-------------------|---------------------------|
|               |                        |                     | (25-500 °C)     | (25-600 °C) |                          |                                       |                                       |                   |                           |
| 1040-1115     | 800                    | 1260                | 14.4            | 14.6        | HV 5<br>b) 205<br>a) 220 | N/mm <sup>2</sup><br>b) 500<br>a) 520 | N/mm <sup>2</sup><br>b) 605<br>a) 630 | %<br>b) 7<br>a) 6 | g/cm <sup>3</sup><br>19.2 |

b: Achievable values after ceramic veneering  
 a: 15 min./950 °C + 15 min./450 °C

\* Measurements performed on samples pursuant to DIN EN ISO 9693

Instructions for Use

Note: Please note the melting range (1040-1115 °C) when setting the firing temperature for dental ceramics and when soldering. Please check and calibrate your ceramic furnace for this temperature range.

Modelling: When modelling and finishing the frameworks, make sure to maintain a minimum wall thickness of 0.3 mm for single crowns and 0.5 mm for abutment crowns.

Sprueing: Fit the waxed-up bridge framework with adequately dimensioned sprues.

Bar casting: The diameter of the head canal is 3–4 mm. The manifold canal has a diameter of 4-5 mm depending on the volume of the casting. The volume of the bar should correspond to that of the casting. The distance between the manifold canal and the waxed-up framework should be 5-8 mm, diameter 3 mm.

Investing: In phosphate-bonded investment compounds (such as Deguvest® SR, F, and CF).

For rapid heating we recommend Deguvest® Impuls and SR.

Wax elimination/preheating: At 300 °C and 800 °C depending on flask size. Follow the manufacturer's recommendations.

Casting: The alloy can be safely processed in resistance-heated furnaces in graphite crucibles, or in open-flame casting units in the ceramic-melting trough.

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Wax elimination/preheating: At 300 °C and 800 °C depending on flask size depending on the manufacturer's recommendations.

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Casting temperature 1260 °C

At least one-third of new material must be used. Sandblast sprues/cones before reuse to thoroughly remove all oxides or residual investment.

Post-heating: Resistance-heated: 60 s, Prestomat: 45 s, oxygen/propane; HF and arc-casting machine: 5-10 s.

Divesting: Allow the flask to cool to room temperature after casting and divest the castings using 50 µm Al<sub>2</sub>O<sub>3</sub> beads at 2 bar of pressure away from the investment.

Do not use a hammer for divesting!

Finishing and oxidation: The material is finished with stagger-toothed tungsten-carbide cutters. Do not use diamond cutters. Sandblast the framework with Al<sub>2</sub>O<sub>3</sub> 110 µm at a pressure of about 2 bar and atmospherically oxidize at 900 °C for 10 min. The oxide is removed by pickling in a Neacid® bath.

Veneering: Dental ceramic materials such as Duceram® Plus and Duceram® Kiss are suitable for ceramic veneering. Observe the ceramic manufacturers' recommendations.

Oxide removal: After ceramic firing or after the final heat treatment is carried out, the oxides should be carefully removed. This is done by pickling with Neacid® or by mechanical removal.

|                |   |                  |           |
|----------------|---|------------------|-----------|
| Soldering:     | Before firing   | Degudent®-Lot G1 | (1030 °C) |
|                | After firing  | Lot DG 750       | (750 °C)  |
| Flux:          | DS 1; T   |                  |           |
| Laser welding: | Laser-welding wire is available in the sizes 0.35 mm. |                  |           |
| Annealing:     | 15 min/450 °C   |                  |           |

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