Degubond® 4

Au and Pt-metals	Au	Pd	Ag	Sn	
78.7	49.6	29.0	17.5	3.0	mass-%

Au and Pt-metals	Ir	Ga	Ta	Re	
78.7	0.1	0.5	0.1	0.2	mass-%

REF
1406 0002

50534415/a

DeguDent

GB

Gold-reduced alloy for the fabrication of dental restorations (for ceramic bonding)

Type extra hard Color white

Contraindications: Do not use in known hypersensitivity to one or several metals contained in the alloy.
Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

Interactions: Avoid occlusal and proximal contacts of different alloy types.

Avoid inhaling dust and vapours while in contact with dental alloys. Use suitable vacuum appliances for protection against dust and vapours. Use a facemask or respiratory protection.

For additional information concerning handling, reuse, dust protection and application notes, please consult our brochure "Precious-metal dental alloys – Processing instructions".

For dental use only.

Melting range	Preheating temperature	Casting temperature	Mean linear CTE		Vickers hardness	0,2% Yield strength*	Tensile strength*	Elongation*	Density
			(25-500 °C)	(25-600 °C)					
1160-1280 °C	850 °C	1400-1450 °C	14,5 µm/mK	14,7 µm/mK	HV 5 b) 230 a) 250	N/mm ² b) 530 a) 550	N/mm ² b) 720 a) 720	% b) 11 a) 9	g/cm ³ 14,5

a: 15 min/950 °C + 15 min/700 °C

b: Achievable values after ceramic veneering

g: Achievable values after slow cooling inside the mould

* Measurements performed on samples pursuant to EN ISO 9693

Instructions for Use

Note: Please note the melting range (1160-1280 °C) when setting the firing temperature for dental ceramics and when soldering. Please check and calibrate your ceramic furnace for this temperature range.

Modelling: When modelling and finishing the frameworks, make sure to maintain a minimum wall thickness of 0.3 mm for single crowns and 0.5 mm for abutment crowns.

Sprueing: Fit the waxed-up bridge framework with adequately dimensioned sprues.

Bar casting: The diameter of the head canal is 3-4 mm. The manifold canal has a diameter of 4-5 mm depending on the volume of the casting. The volume of the bar should correspond to that of the casting. The distance between the manifold canal and the waxed-up framework should be 5-8 mm, diameter 3 mm.

Investing: In phosphate-bonded investment compounds (such as Deguvest® SR, F, and CF).

For rapid heating we recommend Deguvest® Impuls or SR.

Wax elimination/preheating: At 300 °C and 850 °C depending on flask size. Follow the manufacturer's recommendations.

Casting: The alloy can be safely processed in resistance-heated furnaces in ceramic crucibles, or in open-flame casting units in the ceramic-melting trough.

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Casting temperature 1400-1450 °C

At least one-third of new material must be used. Sandblast sprues/cones before reuse to thoroughly remove all oxides or residual investment.

Post-heating: Resistance-heated: 180 s, Prestomat: 165 s, oxygen/propane; HF and arc-casting machine: 5-10 s.

Divesting: Allow the flask to cool to room temperature after casting and divest the castings using 50-µm Al₂O₃ beads at 2 bar of pressure away from the investment.

Do not use a hammer for divesting!

Finishing and oxidation: The material is finished with stagger-toothed tungsten-carbide cutters. Do not use diamond cutters. Sandblast the framework with Al₂O₃ 110 µm at a pressure of about 2 bar and atmospherically oxidize at 980 °C for 10 min. Remove the oxide by once more sandblasting the object with Al₂O₃ 110 µm at a pressure of 2 bar.

Veneering: Dental ceramic materials such as Duceram® Plus and Duceram® Kiss are suitable for ceramic veneering. Observe the ceramic manufacturers' recommendations.

Oxide removal: After ceramic firing or after the final heat treatment is carried out, the oxides should be carefully removed. This is done by pickling with Neacid® or by mechanical removal.

Soldering: Before firing Degudent®-Lot N1W (1070 °C)
After firing Degulor®-Lot 2 (745 °C)

Flux: DS 1, T

Laser welding: Laser-welding wire is available in the sizes 0.35 mm.

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