GOLDEN  
G A T E  
SYSTEM

## Cergo Kiss

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## 1. General notes on Cergo Kiss

### Purpose

Pressable ceramic material for inlays, onlays, veneers and crowns. Cergo Kiss pressable ceramics can be excellently veneered with Duceragold Kiss veneering ceramics.

Cergo Kiss is a component of the GoldenGate System. For dental use only.

### Contraindications

- Fabrication of bridgework
- Deep subgingival preparations
- Bruxism or other types of parafunction

### Technical Data

- Dental ceramic material: Type 2, Class 1 according to DIN EN ISO 6872
- Flexural strength and chemical solubility according to DIN EN ISO 6872

### Precautionary notes for medical products

#### Adverse effects and interactions

We are not aware of any risks or adverse effects related to Cergo Kiss pressable ceramics.

If properly processed and used, adverse effects of these medical products will be highly unlikely. However, reactions of the immune system (such as allergies) or localized paraesthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely excluded in principle. Should you hear or be informed of any adverse effects – even when doubtful – we would like to request notification.

In patient hypersensitivity to Cergo Kiss pressable ceramics or one of its ingredients, this medical product may not be used or only under the particular scrutiny of the dentist or physician in charge.

Known cross-reactions or interactions of this medical product with other medical products or material already present in the oral environment must be taken into consideration by the dentist or physician in charge when selecting this medical product.

Notify the dentist or physician in charge of all factors described above if you use this medical product for a custom construction. When working with these materials, make sure to comply with the Instructions for Use and the pertinent Material Safety Data Sheets (MSDS).

## 2. Processing

### Die preparation (Figure 1)

- Establish the preparation margin.
- Use die hardener on the die.
- Apply die spacer on the dye.
- For veneers, inlays, onlays or crown frameworks.

### Precautions

- Use only glass beads (2–4 bar), not aluminium oxide, for air-abrading the frameworks.
- Use water cooling when working with the pressable ceramics to avoid local overheating.

### Safety notes

- Do not inhale dust particles during grinding.

### Transport and storage

- Liquids: Store containers tightly closed at temperatures above 10°C.
- Protect powder from moisture.

### Symbols on product labels

REF	Product code
LOT	Batch or lot number
	Use before
	Follow Instructions for Use
	Do not reuse
	10°C Lower temperature limit
	Keep dry

### Combinable liquids

- Dentine/incisal, etc.:
  - Ducera® Liquid SD
  - Ducera® Liquid Form
- Paints/glazes:
  - Ducera® Liquid Stain improved

### Adhesion

Follow the manufacturer's recommendations for composite.

Cergo Kiss all-ceramic restorations cannot be temporarily cemented and then worn in situ. Cergo Kiss must be adhesively attached.

In order not to compromise the quality of the definitive cementing, temporaries should only be cemented using non-eugenol cements.

Release date: April 2005

### Wax modelling (Figure 2)

Use only wax materials that burn out without residue. Use Isolit isolating liquid.

For all restorative geometries, the wax modelling must have a minimum thickness of 0.7 mm (maximum wall

thickness: 2.5 mm). For posterior teeth, wax up a reduced anatomical form (including the occlusal profile), a type of dental core, in order to ensure sufficient stability.

#### **Sprueing (Figure 3)**

- Wax up the pressing sprue in the position of maximum circumference, in the direction of the flow (incisal for copings, proximal for inlays). Make sure the diameter of the sprue is not reduced at the contact point; ensure a large connection area.
- Sprue length for Cergo Press: 5–6 mm; for Multimat® 2 Touch+Press: 2–3 mm.
- Position all pressing objects at the edge of the base at the same height, i.e. with all object margins at the same level. This will automatically result in sprues of different length. That is fully acceptable.

#### **Up to 0.6 g wax weight:**

1 pressing pellet 100-g muffle

#### **Up to 1.4 g wax weight:**

2 pressing pellets 200-g muffle

#### **Up to 4.8 g wax weight:**

8 pressing pellets Cercon® PressMaster muffle (400 g)

#### **Investing (Figure 4)**

- Place the muffle ring on the muffle former.
- Mix the investment material (Cergo fit or Cergo fit SPEED) as per the manufacturer's instruction. Vibrate lightly into the muffle avoiding bubble formation, until all objects are completely covered with investment. Top off the muffle without vibrating.
- Place the muffle gauge to ensure correct muffle height and to provide for a 90-degree position of the muffle base relative to the pressing direction of the pressing oven. When placing the muffle gauge, investment compound must exit from the gauge opening.
- Let muffle rest until the investment (Cergo fit or Cergo fit SPEED) has set. Protect from vibration.

#### **Pre-heating (Figure 5)**

- When the investment has set, remove the various parts of the muffle.
- The muffle gauge contact point is cleaned and checked for 90-degree orientation. The muffle must not "rock", as this may result in pressing errors.
- The muffle is placed in the pre-heating oven together with the aluminium oxide pressing die, which must be free of ceramic residue. The base temperature of the pre-heating oven will depend on your choice of investment material (refer to pertinent Instructions for Use).
- Heat the muffle from room temperature to 850 °C (refer to investment manufacturer's recommendations for the proper heating procedure).
- When using Cergo fit SPEED investment, you may place the muffle directly into the oven pre-heated to 850 °C after a setting period of 15 minutes.

#### **Pressing technique**

- Activate pressing program (see Pressing programs on p. 18)
- Remove the muffle from the pre-heating oven once the starting temperature of 700 °C has been reached.
- Place Cergo Kiss pellets in the muffle channel.
- Position the pressing die.
- Immediately place the muffle thus prepared in the pressing oven and start the pressing program (for Cergo press: press the start button one more time).
- Once the pressing process has been completed (after approximately 45 minutes), remove the muffle from the pressing oven immediately and allow to slowly cool to room temperature.

#### **Divesting (Figure 6–7)**

- Indicate the position of the pressed objects using the pressing die.
- Make a deep cut into the investment compound, preferably using a diamond-covered and sintered large carbide disc or a carbide disc for metal castings.
- Separate the part of the muffle containing the aluminium-oxide pressing die from the rest of the muffle by turning in opposite directions.
- Use a jet polisher (50 µm, 4 bar; do not use a revolving-jet unit) to remove the investment all the way to the pressed objects. In doing so, make sure not to touch the objects (abrade the "internal cylinder" and remove the external ring of investment).
- Once the objects have become visible, continue abrading across the area using reduced pressure (2 bar).
- Clean the aluminium-oxide pressing die using aluminium-oxide abrasive and rinse.

#### **Finishing (Figure 8)**

- Cut off the sprues using water cooling and no pressure.
- The restoration to be veneered or painted is then finished using diamond burs or stones. This will maintain the surface texture function.
- Use a steam cleaner (Aquaclean) to clean the object. Do not approach the object too closely with the nozzle (risk of local overheating).

**Cergo Kiss: All pellets are compatible with painting and layering techniques.**

#### **Staining technique (Figure 9)**

- Detailed characterization of the ceramic surface as a stain firing (layering technique).
- Complete colour characterization after pressing with the P1 and P2 pellet. Reproduces all "V" colours.
- Gross or detailed characterization of all pressed pellets (staining and layering technique).
- Use for all CTE-adjusted ceramics.
- In principle, all pellets can be used with body and incisal stains.

### Build-up technique (Figure 10)

For veneers, inlays and crowns, use the layering technique and Duceragold Kiss. Dispense with the opaque firing and start building up the ceramics by analogy with ceramo-metal restorations. Firing a ceramic shoulder onto the Cergo Kiss framework is discouraged for reasons of stability.

#### Please note:

- To ensure the stability of the structure, allow adequate space for the Cergo restoration.
- To avoid tension, eliminate sharp line angles.
- Avoid preparation undercuts.

- The dentist should fill deep cavities by placing a core restoration.
- Pay particular attention to the design of the margin.
- Allow for a uniform wall thickness wherever possible.
- When preparing anterior or posterior teeth, maintain a framework to veneering material ratio of 50 : 50.
- For inlays, onlays and veneers, maintain a framework to veneering material ratio of 75 : 25.

## 3. General firing and pressing recommendations

### Pressing programs

Cergo press	Starting temp. °C	Vacuum	Heating °C/min	Final temp. °C	Holding min	Press tech. min	Pressure bar
100g-/200-g-muffel	700	cont.	60	980	20:00	20:00	4,5
Cercon PressMaster	700	cont.	60	980	40:00	40:00	4,5

ProFire press	Start °C	Heating °C/min	Pressingtemp °C	Holding min	Pressing min	Pressure bar
100g-/200-g-muffel	700	60	980	20:00	20:00	2,7
Cercon PressMaster	700	60	980	40:00	40:00	3,0

Multimat Touch&Press	Starting temp. °C	Vacuum hPa	Heating rate °C/min	Pressing temp. °C	Vacuum time min	Holding time min	Pressing time min	Pressure bar
100g-/200-g-muffel	700	50	60	960	40:00	20:00	20:00	2,7
Cercon PressMaster	700	50	60	960	40:00	40:00	40:00	2,7

### Firing sequence

General firing program	Pre-heating temp. °C	Drying time min	Heating rate °C/min	Firing temp. °C	Holding time min	Vacuum hPa	Long-term cooling
<b>Layering technique</b>							
Dentine 1	450	6:00	55	800	1:00	50	–
Dentine 2	450	5:00	55	800	1:00	50	–
Glaze	450	3:00	55	790	1:00	–	–
Correction	450	5:00	55	740	1:00	50	–
<b>Painting technique</b>							
Paint	450	3:00	55	800	1:00	–	–
Glaze	450	3:00	55	800	1:00	–	–

For your orientation: The lustre of the veneer after the first firing should be comparable with the veneer of precious-metal frameworks of the GoldenGate System (Figure 11: completed crown).

The values listed here are intended for orientation only and should be regarded only as guidelines. Your firing results may differ. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. Therefore, the guideline values will have to be adapted individually for each firing.

We recommend running a test firing cycle to evaluate the performance of the furnace used. We have compiled and checked all values and other data with great care. However, we cannot under any circumstances be liable for your results.

The glaze firing is effected at no more than 800 °C using LFC glaze and body and incisal paints.

For the most current firing recommendations you may also visit [www.kiss-ceramics.com](http://www.kiss-ceramics.com)

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Dentsply LTD  
 DeguDent Division  
 Building 1  
 Aviator Park  
 Station Road  
 Addlestone  
**KT15 2PG**  
**U.K.**  
 Phone: +44 1932 837315

DeguDent  
 Avda de Burgos 39, 2a pl.  
**28036 Madrid**  
**Spain**  
 Phone: +34 900 811 261

DeguDent Moscow  
 Malaya Nikitskaja 27-2, 1st Floor  
**121069 Moscow**  
**Russia**  
 Phone: +7 495 937 5749

DeguDent Ukraine  
 Mechanisatorow str. 2  
**03035 Kiev**  
**Ukraine**  
 Phone: +38 050 4489755

US-Representative  
 Dentsply International Inc.  
 Prosthetics Division  
 570 West College Avenue  
**York, PA 17405-0872**  
**USA**  
 Phone: (800) 877-0020

Dentsply Canada Limited  
 161 Vinyl Court  
**Woodbridge, Ont. L4L 4A3**  
**Canada**  
 Phone: +1 800 263 1437

DeguDent Industria e Comercio Ltda.  
 Rua Funchal, 375-5<sup>th</sup> floor  
 Vila Olimpia  
**São Paulo**  
**Brazil**  
 Phone: +55 11 38 97 6877

Dentsply Australia PTY Ltd.  
 11-21 Gilby Road  
**Mount Waverley Victoria 3149**  
**Australia**  
 Phone: +61 3 95388241

Dentsply Asia  
 25/F, Pacific Link Tower  
 Southmark, 11 Yip Hing Street  
**Wong Chuk Hang**  
**Hong Kong**  
 Phone: +852 28700336

Dentsply Int'l Inc. Taiwan Branch  
 7F-11, 351, Sec. 2, Chung Shan Rd.  
**Chung Ho, Taipei County**  
**Taiwan**  
 Phone: +886 2 2223-7777 Ext.: 11

Dentsply-Sankin KK  
 3 F & 4F, 1-8-10 Azabudai, Minato-ku  
**Tokio 106-0041**  
**Japan**  
 Phone: +81 3 5114 1001

Cergo Kiss	Color samples			Cergo	
<b>C</b>			Transluzenz nimmt zu / increase of translucency	<b>C</b>	T
<b>C0</b>		A1 • B1 • C1		<b>C0</b>	-
<b>L0</b>		A2 • A3 • B2 • C2 • D2 • D3		<b>L0</b>	-
<b>L1</b>		A1 • B1 • B2 • C1		<b>L1</b>	S1
<b>L2</b>		A2 • A3 • C2 • D2 • D3		<b>L2</b>	S2
<b>LY</b>		B		<b>LY</b>	TY
<b>LR</b>		A		<b>LR</b>	TR
<b>LG</b>		C • D		<b>LG</b>	TG
<b>LB</b>				<b>LB</b>	TB
<b>I1</b>		A1 • B1 • C1		<b>I1</b>	T1
<b>I2</b>		A2 • B2 • C2 • D2	<b>I2</b>	T2	
<b>I3</b>		A3 • A3,5 • B3 • C3 • D3	<b>I3</b>	S3	
<b>I4</b>		A4 • B4 • C4 • D4	<b>I4</b>	S4	
<b>P1</b>		A1-D3 / Body 1-8	<b>P1</b>	W1	
<b>P2</b>		D4 / Body 9	<b>P2</b>	W2	
<b>PF</b>		Bleach	<b>PF</b>	F1	

crystal **C** Pellets  
 - Inlays / Veneers  
 - Schichttechnik /  
 Build-up  
 technique

lucent **L** Pellets  
 - Inlays / Veneers  
 - Schichttechnik /  
 Build-up  
 technique


incisal **I** Pellets  
 - Inlays / Veneers  
 - Schichttechnik /  
 Build-up  
 technique

paint **P** Pellets  
 - Schichttechnik /  
 Build-up  
 technique  
 - Maltechnik /  
 Staining  
 technique

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






DeguDent  
A Dentsply International Company

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DeguDent GmbH