# Duceragold® Kiss *improved* Pastenopaker

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# Duceragold® Kiss *improved* paste opaque

- Modified for easier application
- Note changes in firing sequence
- Proven shades are maintained

Duceragold  
Keep it simple and safe.

• NEW • NEW • NEW •

## New opaque firing table for the new and improved Duceragold® Kiss *improved* paste opaque

Duceragold® Kiss/Profire®\*

	Drying		Closing	Preheating		Vacuum on/off			Increase	Firing	Firing	Vacuum	Tempering	Alloy
	°C	min.	min.	°C	min.	Continuous	on	off	°C	°C	min.	min.	°C min.	
<b>PO Opaque 1 + 2</b>	135	6:00	6:00	500	1:00	cont.	500	780	55	780	1:00	1:00		

Multimat® Touch & Press/Multimat® Touch & Press 2/Multimat® NT/NTpress/Multimat®Easy\*

	Preheating	Drying	Preheating	Vacuum	Heating	Firing	Vacuum	Firing	Tempering	Tempering	Cooling
	°C	min.	min.	hPa	°C/min.	°C	min.	min.	°C	min.	Step
<b>PO Opaque 1 + 2</b>	500	12:00	1:00	50	55	780	1:00	2:00			

\* The values listed here are intended for orientation only and should be regarded only as guidelines. Actual firing results may vary. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. The guideline values therefore need to be adapted individually. We recommend a test firing cycle to check the performance of the furnace. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

**DeguDent**  
A Dentsply Company

Duceragold® Kiss/Programat P 90/P 95\*

Duceragold® Kiss – Firing recommendations for Programat P 90/P 95										
	Firing	Standby temperature °C	Closing time min.	Heating rate min.	Firing temperature °C	Holding time min.	Vacuum on °C	Vacuum off °C	Long-term cooling L °C	Alloy
<b>PO Opaque</b>	Paste 1 + 2	500	13:00	55	780	2:00	500	779		

Duceragold® Kiss/Cergo® \*

	Drying °C min.	Closing min.	Preheating °C min.	Vacuum on/off ein aus			Increase °C	Firing °C	Firing min.	Vacuum min.	Tempering °C min.	Alloy
<b>Opaque 1 + 2</b>	135 6:00	6:00	500 1:00	Cont.	500	780	55	780	1:00	1:00		

Duceragold® Kiss/Dekema Austromat 3001\*

Firing	
<b>Opaque 1 + 2</b>	C500 T360 T360 • L9 T60 V9 T055 • C780 T60 V0 T60 C0 L0 T2 C500

Duceragold® Kiss/Dekema D4\*

Program															
Starting temperature °C	Professional														
	Firing	Drying min.	Closing min.	Preheating °C min.	Temperature 1 °C ↑ min.			Temperature 2 °C ↑ min.			Temperature 3 °C ↑ min.			Vacuum °C ↑ min.	
	<b>Opaque 1 + 2</b>	6:00	6:00	500 1:00	780	55	1:00						780	100%	1:00

\* The values listed here are intended for orientation only and should be regarded only as guidelines. Actual firing results may vary. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. The guideline values therefore need to be adapted individually. We recommend a test firing cycle to check the performance of the furnace. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

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**DeguDent GmbH**