Kiss Artist Kit

Ergänzung zur Gebrauchsanweisung
Directions for Use – Supplement
Complément au mode d'emploi
Integrazione delle istruzioni per l'uso
Ampliación de las instrucciones de uso
Aanvulling op de gebruiksaanwijzing

DeguDent
A Dentsply Company

Kiss Artist Kit

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Ergänzung zur Gebrauchsanweisung	4
Directions for Use – Supplement	10
Complément au mode d'emploi	16
Integrazione delle istruzioni per l'uso	22
Ampliación de las instrucciones de uso	28
Aanvulling op de gebruiksaanwijzing	34

Kiss Artist Kit

Action-i Dentine

Corn

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Action-i dentine

Corn

Translucent, similar to opal incisals to support D-shades, also suitable for alternating layers for mamelons

Butter

Translucent, similar to opal incisals to support B-shades, also suitable for alternating layers for mamelons

Honey

Translucent, similar to opal incisals to support A-shades, also suitable for alternating layers for mamelons

Lavender

Translucent, similar to opal incisals to support C-shades, also suitable for alternating layers for mamelons

Creme

Modifier; to be applied in the upper third of the dentine layer. Can be used pure or mixed.

Marble

Modifier; to be applied in the upper third of the dentine layer. Can be used pure or mixed.

Chocolat

Modifier for additional darkening of darker tooth shades; to be applied to the entire body, can be used pure or mixed

Fluorescent dentine

Gray Inside (GI)

A unique fluorescent dentine that can be used to compensate for high-value framework contours in the incisal region

Transpa effects

Transpa Red (TR)

For additional customization of reddish tooth shades

Gums


Gum 3

Complement to Gum 1 and Gum 2 – chromatic

Gum 4

Complement to Gum 1 and Gum 2 – highly chromatic

General firing recommendations

Please observe the instructions of Cercon ceram Kiss, Kiss and Duceragold Duceram Kiss. There you will find information on the intended use, contraindications, warnings, cautions and technical data.

Cercon ceram Kiss

Firing	Preheating °C	Drying time min	Heating rate °C/min	Final temp. °C	Holding time min	Vacuum hPa	Extended cooling min
Paste Liner 1	575	8:00	55	970	1:00	50	–
Paste Liner 2	575	8:00	55	960	1:00	50	–
Margin 1+2	450	6:00	55	850	1:30	50	–
Dentine 1	450	5:00	55	830	1:30	50	–
Dentine 2	450	5:00	55	820	1:30	50	–
Glaze 2	450	3:00	55	800	1:00	–	6:00
Correction (Final Kiss)	450	5:00	55	680	1:00	50	6:00
Final shoulder (F-SM)	450	5:00	55	680	1:00	50	6:00

Slow cooling is mandatory after the last firing; this includes correction firings of restorations after try-in.

The firing temperature must be adapted to the number of units fired in the same cycle. Five to ten units require an increase by 5 °C to 10 °C; more than ten units require an increase by 10 °C to 20 °C.

The values listed here are intended for orientation only and should be regarded only as guidelines. Actual firing results may vary. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. The guideline values therefore need to be adapted individually for each firing. We recommend a test firing cycle to check the performance of the furnace. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

For up-to-date firing recommendations please visit www.Kiss-keramik.de

General firing recommendations

Duceragold Kiss

Firing			Degunorm				
	Pre-heating temp. °C	Drying time min	Heating rate °C/min	Firing temp. °C	Holding time min	Vacuum hPa	Tempering
Oxide firing	575	0:00	55	780	5:00	50	–
Degunorm classic Base	575	7:00	55	780	1:00	50	–
Paste opaque	575	7:00	55	780	1:00	50	–
Shoulder 1	450	5:00	55	780	1:00	50	–
Shoulder 2	450	5:00	55	780	1:00	50	–
Dentine 1	450	6:00	55	780	1:00	50	3 min/720°C
Dentine 2	450	5:00	55	780	1:00	50	–
Glaze 2	450	3:00	55	770	1:00	–	–
Correction (Final Kiss)	450	2:00	55	720	1:00	50	50
Final Shoulder	450	2:00	55	720	1:00	50	–

Tempering or cooling for all firings from the 1st dentine firing on, including glaze firing for [Degunorm eco](#) and [Econolloy Au](#) (temper for 3 min at 720 °C or cool 5 min to 680 °C; alternatively, use recommended cooling level). Without tempering or cooling for [Degunorm logic](#).

The values listed here are intended for orientation only and should be regarded only as guidelines. Actual firing results may vary. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. The guideline values therefore need to be adapted individually for each firing. We recommend a test firing cycle to check the performance of the furnace. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

Note:

Note the following procedure for pre-drying try-in objects made of Duceragold Kiss hydrothermal ceramic material:

- Heat objects to about 80 °C to 90 °C and allow to pre-dry for 30 minutes, or 1 to 1.5 hours for larger objects, until the moisture and organic residue accrued during try-in have dried.
- Depending on the degree of contamination, steam down or sandblast the object lightly with 50 µm aluminium oxide to remove any organic residue.
- Heat the restoration to 450 °C at a heating rate of 5 °C to 10 °C per minute and hold for another hour.
- Conclude by performing any corrections required.

Duceram Kiss

	Firing	Preheating °C	Drying time min	Heating rate °C/min	Final temp. °C	Holding time min	Vacuum hPa	Tempering
	Oxide	For details, please refer to the processing instructions for the respective non-precious alloys.						
Bioalloys	Neutral paste	575	7:00	55	900	3:00	50	–
	Paste opaque	575	7:00	55	900	3:00	50	–
	Powder opaque	575	5:00	55	900	3:00	50	–
Conventional alloy	Paste opaque 1 + 2	575	7:00	55	930	2:00	50	–
	Powder opaque 1 + 2	575	5:00	55	930	2:00	50	–
No extended cooling Examples Degudent Kiss	Shoulder 1	575	7:00	55	920	1:00	50	–
	Shoulder 2	575	7:00	55	920	1:00	50	–
	Dentine 1	575	6:00	55	910	1:00	50	–
	Dentine 2	575	4:00	55	900	1:00	50	–
	Glaze 2	575	3:00	55	890	1:00	–	–
	Add-On	575	4:00	55	880	1:00	50	–
Extended cooling from CTE 14.6 µm /m·K	Final Shoulder	450	4:00	55	660	1:00	50	–
	Dentine 1	575	6:00	55	910	1:00	50	3 min/850 °C
	Dentine 2	575	4:00	55	900	1:00	50	3 min/850 °C
	Glaze 2	575	3:00	55	890	1:00	–	3 min/850 °C

Duceram Kiss – Non-precious alloys

Firing	Preheating °C	Drying time min	Heating rate °C/min	Final temp. °C	Holding time min	Vacuum hPa	Tempering	Extended cooling
Oxide	For details, please refer to the processing instructions for the respective non-precious alloys.							
Bonder	575	7:00	55	980	2:00	50	–	–
Paste opaque	575	7:00	55	950	2:00	50	–	–
Powder opaque	575	5:00	55	950	2:00	50	–	–
Shoulder 1 + 2	575	7:00	55	930	1:00	50	–	–
Dentine 1	575	6:00	55	920	1:00	50	3 min/850 °C	to 600 °C
Dentine 2	575	4:00	55	910	1:00	50	3 min/850 °C	to 600 °C
Glaze 2	575	3:00	55	890	1:00	–	3 min/850 °C	to 600 °C
Add-On	575	4:00	55	880	1:00	50	3 min/850 °C	to 600 °C
(Final Kiss) Final Shoulder	450	4:00	55	660	1:00	50	–	–

Stress relief cooling down to base temperature should be performed independent of the CTE value. Please also follow the alloy manufacturers' recommendations. To compensate for the low thermal conductivity of non-precious alloys, it is recommended to increase the firing temperature or to extend the firing time for dentine firings for all bridges spanning more than five units.

The values listed here are intended for orientation only and should be regarded only as guidelines. Actual firing results may vary. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. The guideline values therefore need to be adapted individually for each firing. We recommend a test firing cycle to check the performance of the furnace. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

Wedge-shaped mould for firing samples

Veneering ceramics are more sensitive to underfiring (physical properties such as CTE and strength are not correctly adjusted) than to over-firing. A distinct sign of overfiring is rounding of the edges, which can be alleviated by slightly reducing the firing temperature (by 5 °C to 10 °C).

Procedure

Create a firing sample from the respective Transpa Clear. For this purpose, transfer the ceramic material mixed with SD Liquid ceramic to the wedge-shaped mould and remove the liquid with a paper towel. Remove excess ceramic material with a spatula and level the surface slightly.

Place the Kiss Transpa Clear firing sample and a thin firing pad placed on the firing pins and position the assembly on the firing tray (Figs. 1 and 2). Alternatively, you may use platinum foil instead of a firing pad.

Fig. 1: Position in the middle of the radius


Fig. 2: Enlarged detail from Figure 1

→ Note

For the reasons mentioned, it is necessary to achieve the correct degree of firing for the veneering ceramic; otherwise, best results cannot be guaranteed, and in exceptional cases (extreme underfiring), cracks may form and flaking may occur. To achieve this objective, the general Kiss firing parameters may need to be adjusted to the specific furnace used (which has its own individual firing characteristics). This is achieved in a simple way by determining the firing intensity using the wedge-shaped mould for firing samples.

! Important

Please use the same firing tray and pins as generally used in the laboratory.

This is the only way to make the result representative of your daily laboratory work (crowns/bridges), as firing trays and pins exert a significant influence on the degree of firing achieved (see Figs. 1 and 2).

Always fire the firing sample wedge at object level. Do not place directly to the honeycomb firing carrier.

Fire the firing sample using the same burning program as specified for the veneer in the Kiss Instructions for Use:

Dentine firing	Preheating °C	Drying time Pre-heating time min	Heating rate °C/min	Final temp. °C	Holding time min	Vacuum hPa	Extended cooling
Duceragold Kiss	450	5:00	55	780	1:00	50	–
Duceram Kiss	575	6:00	55	910	1:00	50	–
Cercon ceram Kiss	450	5:00	55	830	1:30	50	–

! Important

When determining the degree of firing, the furnace must be preheated well.

Assessment

Look at the result.

- **What do the edges look like?**
Is the edge perfectly sharp, or is it rather rounded?
- **What does the lustre look like?**
Is it the perfect intermediate lustre, or is it too high or too low?

Place the firing sample on a checkered piece of paper. If fired correctly, you should see the checker lines through the transparent mass (cf. Fig. 3).


Fig. 3:

Top result: Furnace overfires slightly.

Bottom result: You may continue processing.

Assessment	Procedure
The result of the firing sample test looks good: sharp edges, medium lustre and transparent.	You can start processing Kiss ceramics directly using the firing specifications in the Directions for Use.
The firing sample appears to have a higher lustre and the edges are slightly rounded.	The furnace overfires slightly. Decrease the final temperature somewhat (by 5 °C to 10 °C) and produce a new firing sample (all other firing parameters remain as given in the firing table). Subtract the temperature difference from ALL final temperatures in the firing programs in the Instructions for Use. They now reflect to your individual "correction factor".
The firing sample appears to have a lower lustre.	The furnace underfires. Increase the final temperature by 5 °C to 10 °C and produce a new firing sample (all other firing parameters remain as given in the firing table). Repeat until the firing sample has distinct edges and exhibits medium gloss and transparency. Add the temperature difference to ALL final temperatures in the firing programs in the Instructions for Use (including opaque firings, etc.).

→ Notes

Phenomena such as age-related changes or the manufacturer's calibration standards may affect the temperature of the furnace. This may result in the readout of the furnace display differing from the actual temperature inside the furnace.

This is why the result obtained with the firing sample is so important.

When firing large restorations or a large number of smaller units, more energy is required.

In this case, increase the temperatures of the dentin firings by 10 °C to 20 °C, depending on the size of the objects.

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DeguDent GmbH