Gebrauchsinformation
Instructions for use
Mode d'emploi
Istruzioni d'uso
Instrucciones de trabajo

Trimmer GT

Degussa 

Degussa 

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A. Equipment Supplied

- Trimmer GT (Type 3) with about 1.50 m connector cable
- 1 water connection pressure hose with hose clamps
- 1 water outlet elbow piece
- 1 wrench for grinding wheel

B. Technical Data

Power voltage: Voltage see typeplate
Nominal current: 3.2 A
Power frequency: 50 Hz
Power consumption: 380 W
Speed: 1410 r.p.m.
Grinding wheel: 255 mm diam.
Grinding surface (width x greatest height): 100 x 140 mm
Housing: cast aluminium
Height/width/depth: 305/300/370 mm
Weight: 16.5 kg

C. Preparations for Initial Operation

1. The Trimmer should be set up next to a water outlet so that the elbow piece supplied with it can lead directly into the outlet basin.

Note: if possible, the outlet basin should be fitted with a plaster trap to prevent gradual clogging of the outlet pipes.

2. Using the pressure hose supplied, connect the water inlet fitting on the Trimmer to the shutoff valve on a water line. Secure both connections with hose clamps.
3. Connect the Trimmer to the power line in accordance with VDE regulations relating to electrical installations (or equivalent local regulations). The Trimmer is now ready for operation.

D. Operation

1. The Trimmer is equipped with a solenoid valve for automatically opening the feed line for the water used to wash the grinding wheel. The amount of water used for this purpose is adjusted with the shutoff valve at the water line. In order to do this, switch on the Trimmer, the grinding wheel starts to turn, and the solenoid valve opens the water feed line inside the unit. Now open the shutoff valve at the water line until the rotating grinding wheel shows a wet surface. Too much water directed onto the grinding wheel will result in undesirable splashing, and the actual washing effect will not be improved.
2. The support table in front of the grinding wheel can be continuously adjusted, so that model sockets with a conical taper can be ground. The support table is positioned at a right angle with respect to the grinding wheel when the marking on the guide plate is adjusted to coincide with the upper edge of the locking plate.
3. The angular milled marks on the support table are provided to permit symmetrical grinding of the models.
4. The relatively large grinding surface also permits the simultaneous grinding of 2 models in occlusion position.
5. When the grinding process has been completed, the Trimmer should be allowed to run for about 10–15 seconds longer with the water continuing to flow so that all grinding residues are washed out from inside the instrument.
6. When the Trimmer is switched off, the solenoid valve automatically blocks the flow of additional water. The quantity of water set at the shutoff valve on the water line remains unchanged as long as the shutoff valve is not closed after each grinding process. It is advisable to close this shutoff valve overnight.

E. Maintenance

The grinding wheel in the Trimmer can be used on both sides. In order to turn this wheel around or to replace it, first disconnect the power cable plug, then unscrew the 4 star grip screws on the back of the grinding wheel housing, and remove the front section of this housing. Release the hexagon nut with the wrench supplied, and turn or replace the wheel. Reassemble the unit by following these steps in the opposite sequence, being sure that the housing seal is well seated.

If deposits of plaster are found on the part of the grinding wheel housing that has been removed, it is advisable to take advantage of this opportunity to wash this part of the housing thoroughly under running water.

The motor and bearings in the Trimmer do not require maintenance.

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