

Cercon ht – Fully contoured safety and so much more.




cercon[®]
smartceramics

Cercon ht – The zirconia concept for the dental laboratory and the dental surgery.

DeguDent
A Dentsply Company


Cercon ht

Safe, conservative preparation and...



Cercon ht – outstanding translucency and conservative preparation option for fully anatomic crowns and bridges.

Safety of fully anatomic designs under real-life conditions



All previous studies on the wear behaviour of fully anatomic crowns and antagonists have been performed on a glaze or polish base and do not take into account any subsequent occlusal adjustments. Now, for the first time, a scientific study has examined adjusted and repolished Cercon ht crowns. Cercon ht crowns, designed to fully anatomic contour, are demonstrably safe for their antagonists.

Source: University of Regensburg, Germany, 2011

Product names marked ® are registered trademarks of the respective trademark owner(s).

... for multiple indications



Cercon ht – whether you go for fully contoured and conservative or for partially or fully veneered, using the pre-coloured Cercon ht light or medium – the new Cercon stains and the special Cercon ht dentines for Cercon Kiss and Cercon love perfectly highlight the excellent translucency and aesthetics of your high-end restorations.



Multiple indications for Cercon ht



High-end aesthetics with love ht dentines



Build-up with chromatic Cercon ceram Kiss ht dentines



Fully contoured using the staining technique



Combination of several veneer types

Cercon ht

Conservative in the dental surgery

Crown

Conservatively prepared, fully contoured crowns




Cercon ht crowns – a harmonious appearance




Cercon ht crowns 25 – 27

Conservative preparation fully contoured crowns



Classic preparation fully contoured crowns

Minimum wall thickness 1.5 – 2.0 mm



Bridges

Conservatively prepared, fully contoured posterior bridges




Cercon ht, bridge 45 – 47




Cercon ht, bridge 15 – 17

Conservative preparation fully contoured bridges




Classic preparation fully contoured bridges



Cercon ht

The zirconia concept for the dental laboratory


Examples of the many different indications and performance levels in all-ceramic prosthodontics using Cercon ht:






FC

fully contoured



PV

partially veneered



FV

fully veneered



bridge

PV


PV

PV

contoured,
bite



Primary crown



Custom abutment
for titanium base



Custom
abutment with
crown for
titanium base

xilla

Fascination Prosthetics

28275/1202/IZ
Last revision: 02/2012

