



Simply Faster  
to the Finish

Cergo fit Speed –  
Special Investment for  
Pressable Ceramics

**DeguDent**


*A Dentsply International Company*

# Save Time, Win with Good Results: The Rapid Heating Investment Cergo fit Speed.

Cergo® fit Speed is the investment of the next generation. That is because this new type of phosphate-bound graphite- and plaster-free investment reduces the firing process by around 3 hours.

Cergo® fit Speed is the continuous development of our product range and opens new paths to reduce costs by enormous time and energy savings.

Comparison Cergo fit Speed - Conventional Investment  
Time Saving: 3 hours



Cergo fit Speed has everything you need for the best pressing results:

## Cergo® fit Speed – Qualities which you can see.



### Forms of Delivery

Box with 65 x 100 g bags

1,450 ml bottle with Cergo fit Speed Liquid

### Technical Data

Mixing Ratio Powder/Liquid

Setting Time (acc. to Vicat)

Processing time span

Setting Expansion (100 % Liquid)

Thermal Expansion (100 % Liquid)

Compressive Strength

Tests according DIN EN ISO 7490

- Safe processing from mixing to sandblasting
- Creamy flowability by coordinated materials
- No limits regarding shock heating
- Up to three hours time saving during preheating
- High temperature resistance
- Optimal fitting results
- Precise reproduction of details
- Excellent surfaces
- Reproducible results

### REF

53 2532 0200

53 2540 0117

100 g: 25 ml

8–11 min

5–7 min

approx. > 2.0 %

> 1.0 %

5–7 MPa