


MultimatNT

MultimatNTpress



# Multi-talented

Freedom is the total of all possibilities.

**DENTSPLY**

## A multitalented individualist

To improve something that is already good means running a risk. DeguDent is courageous enough to take this challenge when redesigning its Multimat series – since the models of the previous series were widely considered innovative masterpieces, and therefore became the world's most sold ceramic furnaces. For the new series, the developers and dental technicians of DeguDent have made sure to retain existing highlights of excellence while adding useful and future-proof features.

There is one thing we were not able to develop into the MultimatNT – YOUR own personal touch. You are the only one who knows what suits your working style best. You will be impressed of the number of possibilities for customization – configure your MultimatNT so it behaves exactly as you need it. Its intelligent menu-driven interface allows your work with the new MultimatNT a simple, convenient and safe experience.

### MultimatNT

For firing

### MultimatNTpress

For firing and for pressing

- Precise firing thanks to state-of-the-art control technology
- High safety thanks to permanently monitored measuring and control loops
- Compact design – small footprint
- Firing curve is shown graphically on the touch screen and can be configured individually
- Self-explanatory menu and touch screen interface
- Interface icons for immediate access to common firing schedules
- Gigantic memory for 1000 stored firing schedules, allowing ceramic types to be stored in sets of 100
- Data transfer to and from PC by USB stick
- Display of images and video clips on a 5.7-inch TFT screen
- Saving, managing and printing of firing tables and firing logs, simply and easily via a PC
- Maintenance-friendly
- Nearly unlimited data retention after power-out



MultimatNT & MultimatNTpress

The high-resolution 5.7 inch touch screen with the intuitive menu-driven interface.

Excellent connectivity via a PC/USB interface, for saving and managing data etc.

Transfer your data via a USB stick – an up-to-date and future-proof technology.

The system's graphic capabilities ensure a brilliant display of all patient-specific and firing-cycle data.

Technical specifications may change without notice  
Estimated Product Launch 4th Quarter 2007

**DeguDent GmbH**

A Dentsply International Company  
P.O. Box 1364 · 63403 Hanau  
GERMANY  
[www.dentsply-degudent.com](http://www.dentsply-degudent.com)

**DENTSPLY**