



Multivac eco

economical and effective

DeguDent
A Dentsply International Company



The Multivac eco offers everything you need for mixing in everyday laboratory life: a powerful mixer and an efficient vacuum pump.

The economical price/performance ratio and the solid construction make this vacuum mixer the most cost-effective machine in your laboratory. Motor and pump are, of course, maintenance-free and impress with their longevity and everyday usability.

Due to the simple, easy-to-use operation, investment materials and plasters are mixed homogeneously and processed safely.

Due to our technique of evacuating directly through the cover of the mixing bowl, an additional suction hose is unnecessary. All the latest DeguDent mixing bowls are compatible (except the extra large mixing bowls).

DeguDent GmbH
P.O. Box 13 64 · 63403 Hanau
GERMANY
www.degudent.com

Technical Data

Power supply:	220 – 240 V, 50 Hz
Power Consumption:	150 W
Rev:	330 1/min
Required fuse protection:	2 AT
Height x width x depth:	505 x 390 x 315 mm
Weight:	8.8 kg

Supply Form

REF

Vacuum Mixer (wall-mounted) with	
1 Power cord	
1 Replacement Filter	
1 Drilling template	
1 Kitt for wall-mounting	5316 0001

Accessories

REF

Mixing Bowl extra small (50 cm ³)*	5314 0113
Mixing Bowl small (150 cm ³)*	5314 0100
Mixing Bowl large (420 cm ³)*	5314 0101
Replacement Filter	5314 0105

* DeguDent Mixing Bowls are delivered together with cover and stirring mechanism.

DeguDent
A Dentsply International Company