


Application reference

The Solaris® gold bath can be used together with the Solaris® activator in a Heraeus CL-GF unit.


Care must be taken to ensure the mixing ratio of 100 ml of gold bath to 5 ml of activator is maintained.

The amount of the gold bath required for precipitation can be determined by consulting Heraeus' table for the Heraeus CL-GF.

All other device parameter settings can be determined by consulting the original instruction manual of the unit.



Dr. Udo Schusser
Director, Research and Development



Dr. Markus Vollmann
Director, Application Technology